


hawo




## hd 650 D/DE EcoPak

### OPERATING INSTRUCTIONS

Valid for device versions hd 650 D/DE EcoPak

<b>1</b>	<b>INTRODUCTION</b>	<b>3</b>
1.1	PREFACE	3
1.2	LEGEND	3
1.3	IMPORTANT NOTICE	4
1.4	CLEANING	4
1.5	SAFETY INSTRUCTIONS	5
<b>2</b>	<b>BEFORE STARTING</b>	<b>7</b>
2.1	INTENDED USE	7
2.2	DESIGN AND FUNCTIONS	8
2.2.1	Sealing process sequence	8
2.3	INSTALLATION	9
<b>3</b>	<b>BASIC FUNCTIONS</b>	<b>10</b>
3.1	OPERATION	10
3.2	OPERATION	11
3.3	PROCESS VARIABLES	11
3.4	SETTING THE TEMPERATURE CONTROLLER	12
<b>4</b>	<b>TROUBLESHOOTING AND MAINTENANCE</b>	<b>13</b>
4.1	TROUBLESHOOTING CHECKLIST	13
4.2	OFFSET – SETTING THE TEMPERATURE	14
4.2.1	Revoke key lock	14
4.2.2	Offset – Make setting	14
4.2.3	Activate key lock	14
4.3	HAWO CUSTOMER SERVICE	15
4.4	SERVICING PLAN	15
4.5	SPARE PARTS SERVICE	15
4.6	SPARE PARTS ORDERING – ALLOCATION OF ARTICLE NUMBERS	17
4.6.1	Upper sealing die 1.616.024	17
4.6.2	Lower sealing die 1.616.025	17
4.6.3	Upper guiding die 1.619.014	18
4.6.4	Lower guiding die 1.619.015	18
4.6.5	Comprehensive overview	19
<b>5</b>	<b>TECHNICAL DATA</b>	<b>20</b>
5.1	SWITCHING AND WIRING DIAGRAM 230V~	20
5.2	SWITCHING AND WIRING DIAGRAM 115V~	21
5.3	SPECIFICATIONS	22
<b>6</b>	<b>DECLARATION OF CONFORMITY</b>	<b>23</b>



# 1 Introduction

## 1.1 Preface

First of all we would like to thank you for purchasing this sealing device. In these instructions you will find information about using the device, servicing and care as well as process validation.

The sealing device is a microprocessor controlled rotary sealer with a printer for packaging sealable transparent pouches and reels (med. packaging).



Please read these operating instructions carefully before commissioning so that you are familiar with the capabilities of the device and you can make optimum use of its functions.

## 1.2 Legend

	The exclamation mark in the triangle draws your attention to important notes in the operating instructions, which must absolutely be observed.
	This warning sign refers to measures that could result in danger to human health if they are not observed. It is compulsory to observe it.
	Tips with a hand symbol next to them relate to daily practice.



### 1.3 Important notice



In accordance with the intended use, the CE marking is displayed based on the following EU directives: 2006/42/EC, 2014/30/EU and 2011/65/EU.

The Medical Device Directive 93/42/EEC and Medical Device Regulation 2017/745 (MDR) is not applicable to sealing devices.

The limit values of IEC 60601-1 must not be applied in repeated electrical inspections.

The manufacturer accepts no liability whatsoever for damage caused by tests in accordance with standards not listed in the declaration of conformity.

In the event of conversion work or interventions to the device undertaken without the express written permission of the manufacturer, the warranty shall be deemed void and any liability for physical or material damage shall be transferred to the operator.

**Note**

We are constantly improving our products, therefore we reserve the right to modify these operating instructions and the functions described in them.

**These operating instructions apply to products from the hd 650 D/DE EcoPak series.**

### 1.4 Cleaning

Before cleaning, disconnect the mains plug from the socket and disconnect the device from the power supply with the plug.

Clean the device only with a dry or damp soft cloth and a mild cleaning agent. (E.g.: isopropanol, spirit, etc.) Do not allow any water to find its way into the device.

**Caution! Never wet clean the device!**

## 1.5 Safety instructions



1. Our products are in a flawless condition in terms of safety technology when they leave the plant.
2. To maintain this condition, the contents of these safety instructions as well as type plates, labelling and safety instructions attached to the device must be observed while handling the device (transport, storage, installation, commissioning, operation and maintenance).
3. This device is suitable for processing laminated films in the heat-sealing process. See also chapter 2.1 "Intended use".
4. Please check the packaging, and lodge a complaint for any damage with the carrier or parcel service immediately, before installing the device.
5. Before commissioning, ensure that the device does not show any evidence of damage. In case of doubt, contact the manufacturer or a service partner authorised by the manufacturer.
6. Do not operate the device if the power cable or the power plug is damaged. Do not use the device if it does not operate correctly or it is damaged in any way. If the mains cable or the device have been damaged, the device must be repaired by the manufacturer or by one of the manufacturer's authorised service partners.
7. The device must be connected using the mains cable included in the scope of delivery to a protective contact socket with a stable voltage. Operation on IT networks is not permitted.
8. Place the device on a stable base.
9. The device must not be installed or operated in potentially explosive areas.
10. If the sealing device is brought directly from a cold environment into a warm environment, condensation may form. Wait until temperature equalisation has taken place.  
**Starting up the device when it contains condensation causes danger to life!**
11. Repairs and the replacement of wear parts / spare parts must be performed only by the manufacturer or by one of the manufacturer's authorised service partners.
12. Switch off the device when it is not in use, or remove the power plug from the socket.
13. **Before cleaning: Disconnect from the mains! Clean the device only with a dry or damp soft cloth and a mild cleaning agent. Do not allow any water to find its way into the device. Caution! Never wet clean the device!**
14. Do not insert pointed or flat items into the import slot of the device. This can result in damage to the device and instruments.
15. Do not insert items into the louvres of the device. You may receive an electric shock or the device could be damaged.
16. Do not use the device if you have any doubts about its safety.
17. The device must not be installed or operated by persons under 16 years of age.

18. The device must not be operated unsupervised.
19. It is forbidden to operate the device under the influence of drugs or alcohol.
20. Keep hair, clothing and gloves away from moving parts! Loose clothing, jewellery or long hair can be caught by moving parts.



21. The device contains valuable materials that can be recycled and reused. The device should therefore be disposed of at a public disposal facility near you. The device has been labelled in accordance with Directive 2002/96/EC (WEEE) on waste electrical and electronic equipment. This directive governs the return and recycling of scrap equipment within the EU.

HENDERSON  
BIOMEDICAL  
www.henderson-biomedical.co.uk

## 2 Before starting

### 2.1 Intended use

#### SEALING MATERIALS


Sealable paper pouches in accordance with EN ISO 11607-1/EN 868-4	x
Sealable pouches and tubes in accordance with EN ISO 11607-1/EN 868-5 made of film and paper as per EN 868-3	x
Sealable pouches and tubes in accordance with ISO EN 11607-1/EN 868-5 made of film and uncoated materials made of polyolefins as per EN 868-9 (e.g. Tyvek <sup>®1</sup> )	x <sup>2</sup>
Sealable pouches and tubes in accordance with ISO 11606-1/EN 868-5 made of PP fleece or PP non-woven	x <sup>2</sup>

#### NON-SEALABLE MATERIALS


Soft PVC films
Polyamide films
Coated HDPE
Aluminium-laminate film
Polyethylene films
Hard PVC films
Polypropylene films

<sup>1</sup> Tyvek<sup>®</sup> is a registered trademark of E.I. du Pont Nemours.

<sup>2</sup> Release and / or test required



## 2.2 Design and functions



### 2.2.1 Sealing process sequence

- Step 1:** After the sterilisation packaging has been inserted, the feed is automatically switched on.
- Step 2:** The sterilisation packaging is now transported and the area is heated to the set sealing temperature.
- Step 3:** The sealing seam, which is now heated, is pressed together by the sealing rollers and sealed.
- Step 4:** The finished sterilization packaging is transported to the extraction side.
- Step 5:** If no item to be sealed is fed in, the feed switches off after approximately 30 seconds.



2.3 Installation



The device must not be installed or operated in potentially explosive areas.


Only use sockets that are equipped with a protective conductor and where the mains voltage is stable.



The device may only be installed in a dry environment. Heavy dust, steam, dripping water or splashing water impair the function of the device.


Please make sure that the operating voltage corresponds to the specifications on the device's type plate.

- ❶ Please do not transport the device on the peel edge setting.
- ❷ The distance from the device to a wall must be at least 200 mm!



### 3 Basic functions

#### 3.1 Operation




Schritt 1: Plug the mains cable into the mains connection.

Schritt 2: Switch on the device with the mains switch in position "1".  
The indicator light in the switch is lit.

Schritt 3: Set the required sealing temperature on the temperature controller as described in Chapter 3.3.

Schritt 4: The device is heated up and ready for operation as soon as the set sealing temperature is displayed.

## 3.2 Operation



- Schritt 1: Set the required sealing edge width. After loosening the knurled screw, the lower insertion plate can be continuously adjusted for sealing edge widths of 0 - 30mm.
- Schritt 2: Insert packaging into the device from the left via the insertion plate.
- Schritt 3: Remove the sealed package on the outlet side and allow to cool down briefly.



### Check of the sealing seam

If leaks appear, the sealing temperature must be increased. If the film melts, the set temperature is too high.

As per DIN 58953 -7, the suitable sealing temperature is to be determined by test sealing.


## 3.3 Process variables

### Sealing temperature

The temperature is monitored electronically by means of a temperature sensor.

If this varies from the set value by 5°C (requirements set out in DIN 58953-7), the drive is locked.

### 3.4 Setting the temperature controller



#### Input of the set temperature

The 5 key is pressed once

Confirm "MENU" [User Level] selection level by pressing the 5 key twice

Setting the temperature with the 2 [+1°C]      3 [-1°C] keys

The input is activated with the 5 key      → OK appears for 2 seconds

Exit the input by twice pressing the key 4

#### Default setting:

Set temperature      180 °C

## 4 Troubleshooting and maintenance

### 4.1 Troubleshooting checklist

Malfunction	Possible cause	Remedy
Device does not switch on	Power supply	
	Power cable not plugged in	Check the mains connection and, if necessary, plug into a different power outlet.
	Power cable defective	Replace power cable
	Line fuse	Replace line fuse. If the fuse blows again, it is imperative to have the device tested!
Device fails to heat up	Set temperature is too low	Increase set temperature
	Overtemperature fuse has tripped	Push in the overtemperature fuse. If it still trips it is imperative to have the device tested!
	Temperature controller	Replace temperature controller
	Heating cartridge	Check heating cartridges and replace if necessary
	Solid state relay	Replace solid state relay
No material transport	Set temperature not reached	
	Motor faulty	Replace motor
Uneven material transport	Microswitch faulty	Replace microswitch
	Toothed belt damaged	Replace toothed belt
	Toothed belt does not transport	Check the tension of the toothed belt
Loud running noise	Motor faulty	Replace motor
Sealing seam will not hold	Sealing temperature too low	Increase sealing temperature
	Distance between the sealing dies is too large	Set the distance between the sealing dies to 0.5mm
Sealing seam distorted or melted	Sealing temperature too high	Reduce sealing temperature

## 4.2 Offset – setting the temperature

### 4.2.1 Revoke key lock



1. Switch on the device
2. Check if the device is heating
3. Press and hold the two keys "Arrow down" & "Menu" simultaneously.
4. Change level from CONFIG. and PARAM. to NONE.
5. Press the Menu/OK key
6. Controller unlocked

### 4.2.2 Offset – Make setting



1. Press the Menu / OK **1** key
2. MENU level  
Switch to **CONFIGURATION** with the arrow keys
3. Press the Menu / OK **1** key
4. CONFIGURATION level  
Switch to **ANALOG INPUT** with the arrow keys
5. Press the Menu / OK **1** key
6. ANALOG INPUT level  
Switch to **MEASURING VALUE OFFSET** with the arrow keys
7. Press the Menu / OK **1** key
8. Enter offset value  
E.g.: TSET: 180°C  
TACTUAL: (measured) 184°C  
MEASURING VALUE OFFSET input: **+4°C**
9. Press the Menu/OK **1** key
10. Press the Return key twice to return to the display mode.
11. Check values and repeat point 4.2.2 if necessary.

### 4.2.3 Activate key lock



1. Switch on the device
2. Check if the device heats up
3. Press and hold the two keys "Arrow down" & "Menu" simultaneously.
4. Change level from NONE to CONFIG. and PARAM.
5. Press the Menu/OK key
6. Controller locked

4.3 hawo customer service



Your hawo customer service is available from Mon-Fri 07:30 - 16:30 CET at the following number: +49 (0)6261-9770-0.

4.4 Servicing plan



Like all technical devices, your device is subject to technical wear. In order to ensure maximum availability, your device should be maintained by a trained professional at least once a year.

Maintenance cycle	Cleaning	PTFE tape Guiding die	PTFE tape Sealing die	Pressure roller	Toothed belt	Distance of Sealing die
At least every 3 months						
Depending on usage, at least once a year						

Legend:



Check



Replace



Adjust

4.5 Spare parts service



Convenient spare parts order by fax!  
To do this, simply copy the order form.  
The form is on the next page.

- Enter address, fax number and order number
- Enter serial number
- Enter the device type
- Mark required items
- Enter required quantity
- Sign order
- Fax order

→ 444.444  
→ hd 650 DE

Bitte diese Daten bei Ersatzteilbestellung angeben.  
Please state this data when ordering spares.  
Veuillez indiquer ces données en cas de commande de pièces de rechange.

**hawo** 2 990 009

**Gerätebau GmbH**  
D-74847 Obriegheim  
Tel. 06261-620 16  
Fax. 06261-620 15  
Made in Germany

**To:**  
**hawo GmbH**  
**Fax no. +49 (0)6261 977069**


**Sender:**  
 \_\_\_\_\_  
 \_\_\_\_\_

<b>Your order no.</b> _____		<b>Date</b> _____	
<b>Device type</b> _____		<b>Serial number</b> _____	
<input checked="" type="checkbox"/>	<b>Designation</b>	<b>Art. no.</b>	<b>Qty.</b>
<input type="checkbox"/>	Pressure roller	2.230.026	
<input type="checkbox"/>	Gear motor           230 V 115 V	1.212.005 1.212.014	
<input type="checkbox"/>	Time controller       230 V 115 V	1.540.056 1.540.057	
<input type="checkbox"/>	Heating cartridge   115 V/200W	6.536.032	
<input type="checkbox"/>	Fan                     230 V 115 V	6.212.019 6.212.021	
<input type="checkbox"/>	PTFE tape heating die	6.105.125	
<input type="checkbox"/>	PTFE tape guiding die	6.105.138	
<input type="checkbox"/>	IEC cable supply with switch	6.562.009	
<input type="checkbox"/>	Safety element for 6.562.009	6.562.049	
<input type="checkbox"/>	Mains cable 230V	6.593.013	
<input type="checkbox"/>	Mains cable 115V	6.593.014	
<input type="checkbox"/>	Opto - Sensor Motor Start   (650 DE only)	1.561.016	
<input type="checkbox"/>	SST relay	6.460.001	
<input type="checkbox"/>	Temperature sensor	6.564.040	
<input type="checkbox"/>	Digital temperature controller 110 -240 V	6.564.050	
<input type="checkbox"/>	Overtemperature protection	6.564.018	
<input type="checkbox"/>	Toothed belt transport	6.271.001	
<input type="checkbox"/>	Motor – Start switch   (650 D only)	6.562.009	
<input type="checkbox"/>			




4.6 Spare parts ordering – allocation of article numbers


4.6.1 Upper sealing die **1.616.024**




4.6.2 Lower sealing die **1.616.025**



4.6.3 Upper guiding die **1.619.014**




4.6.4 Lower guiding die **1.619.015**




HENDERSON  
BIOMEDICAL  
www.henderson-biomedical.co.uk

4.6.5 Comprehensive overview




## 5 Technical data

### 5.1 Switching and wiring diagram 230V~



0	Temperature controller	6.564.050
1	Temperature sensor	6.564.040
2	Overtemperature protection	6.564.018
3	SST relay	6.460.001
4	115V/200W heating cartridges	6.536.024
5	Gear motor	230V 1.212.005 115V 1.212.014
6	Optical sensor	1.561.016
7	Time relay	230V 1.540.056 115V 1.540.057
8	IEC cable supply with switch	6.562.009
9	Fan	230V~ 6.212.019 115V 6.212.021

5.2 Switching and wiring diagram 115V~



HENDERSON BIOMEDICAL

www.henderson-biomedical.co.uk

5.3 Specifications

**Connection data**

Mains connection	[ V ]	115 / 230
Mains frequency	[ Hz ]	50 / 60
Power consumption max.	[ W ]	390
Mains fuse 115V / 230V	[ A ]	6.3A / 4A

**Mechanical system**

Dimensions	Length	[ mm ]	505
Including infeed section	Width		250
	Height		145
Housing			Metal, powder-coated
Weight	[ kg ]		15
Seal distance from the edge	[ mm ]		0 – 35
Seal seam width	[ mm ]		12
Sealing system			Multi – Line
Length of sealing seam	[ mm ]		Unlimited
Distance from medical product	[ mm ]		>30 (as per DIN 58953-7)

**Process variables/sealing parameters**

Sealing temperature max.	[ °C ]	220
Sealing temperature switch-off tolerance	[ °C ]	±5
Throughput speed [fixed]	[ m / min ]	10
Temperature control tolerance	[ % ]	±2



**Electronics and communication systems**

System	electronic
Electrical protection class	1

**Environmental parameters**

Ambient temperature	[ °C ]	5-25
Heat emission	[ kJ/s ]	0.1
Noise intensity as per Machinery Directive 2006/42/EC Appendix I 1.7.4.2 u.)	[ dB/ A ]	<70
Relative humidity	[ % ]	30-80 non-condensing

## 6 Declaration of conformity

 74847 Obrigheim / Germany	Konformitätserklärung – Declaration of Conformity Déclaration "CE" de Conformité Declaración de conformidad de la C.E. Dichiarazione di conformità - Declaração de conformidade	<b>9.693.004C</b>																				
Gültig ab: 06.08.2020 Valid from:		Seite 1/1 Version 3.05																				
<p>Hiermit erklären wir, daß die Folienschweissmaschinen:                  Herewith we declare that the Foil sealing unit:                  Par la présente, nous déclarons que la gamme de Soudeuse de films plastique:                  Por la presente certificamos que las máquinas embolsadoras modelos:                  Dichiariamo con la presente che le macchine per saldatura di fogli:                  Por este meio se declara que as máquinas de selagem de folhas de plástico:</p> <p style="text-align: center;"><b>hd 650 D ECOPAK / hd 650 DE ECOPAK</b></p> <p>folgenden einschlägigen Bestimmungen und harmonisierten Normen entsprechen:                  complies with the requirements of the following regulations and harmonised standards:                  corresponde aux dispositions suivantes et standards harmonisés:                  objeto de esta Declaración cumple con las siguientes disposiciones:                  Sono conformi alle seguenti disposizioni in materia nonché alle seguenti norme armonizzate:                  correspondem às seguintes determinações e normas harmonizadas:</p> <table border="0"> <tr> <td>EG - Maschinenrichtlinie Machinery directive Directive "CE" rel. aux machines Directiva de Maquinaria de la CE Direttiva CE sulle macchine nella versione Directiva da UE relativa a maquinaria</td> <td>2006/42/EG</td> <td></td> <td></td> </tr> <tr> <td>EMV-Richtlinie Directive CEM Direttiva CEM</td> <td>EMC-directive Directiva de CEM Directiva CEM</td> <td>2014/30/EU</td> <td></td> </tr> <tr> <td>WEEE-Richtlinie Directive WEEE Direttiva WEEE</td> <td>WEEE-directive Directiva de WEEE Directiva WEEE</td> <td>2012/19/EU</td> <td></td> </tr> <tr> <td>RoHS-Richtlinie Directive RoHS Direttiva RoHS</td> <td>RoHS-directive Directiva de RoHS Directiva RoHS</td> <td>2011/65/EG</td> <td></td> </tr> <tr> <td>Harmonisierte Normen Standard harmonise Norme armonizzate</td> <td>Harmonized standards Las normas armonizadas Normas harmonizadas</td> <td>EN ISO 12100:2010 EN ISO 13857:2008</td> <td>EN 60204-1:2018 EN 61000-6-1:2019 EN 61000-6-3/A1:2012</td> </tr> </table> <p>Verantwortliche Person für die Technischen Unterlagen siehe unten                  Responsible person for technical documentation see below                  La personne responsable pour la documentation technique est mentionnée au-dessous</p> <div style="text-align: center;">                   hAWO GmbH                  74847 Obrigheim / Germany                  Torsten Ehrhardt                  +49 (0) 6261 / 9770-0                  info@hAWO.com                  Prokurist / authorized officer             </div> <p>hAWO GmbH, Obere Au 2, D-74847 Obrigheim, Germany</p>			EG - Maschinenrichtlinie Machinery directive Directive "CE" rel. aux machines Directiva de Maquinaria de la CE Direttiva CE sulle macchine nella versione Directiva da UE relativa a maquinaria	2006/42/EG			EMV-Richtlinie Directive CEM Direttiva CEM	EMC-directive Directiva de CEM Directiva CEM	2014/30/EU		WEEE-Richtlinie Directive WEEE Direttiva WEEE	WEEE-directive Directiva de WEEE Directiva WEEE	2012/19/EU		RoHS-Richtlinie Directive RoHS Direttiva RoHS	RoHS-directive Directiva de RoHS Directiva RoHS	2011/65/EG		Harmonisierte Normen Standard harmonise Norme armonizzate	Harmonized standards Las normas armonizadas Normas harmonizadas	EN ISO 12100:2010 EN ISO 13857:2008	EN 60204-1:2018 EN 61000-6-1:2019 EN 61000-6-3/A1:2012
EG - Maschinenrichtlinie Machinery directive Directive "CE" rel. aux machines Directiva de Maquinaria de la CE Direttiva CE sulle macchine nella versione Directiva da UE relativa a maquinaria	2006/42/EG																					
EMV-Richtlinie Directive CEM Direttiva CEM	EMC-directive Directiva de CEM Directiva CEM	2014/30/EU																				
WEEE-Richtlinie Directive WEEE Direttiva WEEE	WEEE-directive Directiva de WEEE Directiva WEEE	2012/19/EU																				
RoHS-Richtlinie Directive RoHS Direttiva RoHS	RoHS-directive Directiva de RoHS Directiva RoHS	2011/65/EG																				
Harmonisierte Normen Standard harmonise Norme armonizzate	Harmonized standards Las normas armonizadas Normas harmonizadas	EN ISO 12100:2010 EN ISO 13857:2008	EN 60204-1:2018 EN 61000-6-1:2019 EN 61000-6-3/A1:2012																			
hAWO GmbH Obere Au 2-4 74847 Obrigheim / Germany	T + 49 (0) 6261 / 9770-0 F + 49 (0) 6261 / 9770-69 info@hAWO.com www.hAWO.com	Amtsgericht Mannheim: HRB 441011 Geschäftsführer: Christian Wolf Firmensitz: Obrigheim	This document and the contents hereof are considered proprietary and confidential information of hAWO and disclosure to unauthorized individuals or dissemination, publication, or copying is prohibited without prior written consent by hAWO GmbH, 74847 Obrigheim, Germany.																			

9.693.004 Version 2.01



**hawo**



**hawo GmbH**  
Obere Au 2–4  
74847 Obrigheim  
Germany  
T +49 (0) 6261 / 9770 - 0  
F +49 (0) 6261 / 9770 - 69  
info@hawo.com  
www.hawo.com

**hawo USA**  
321 North Clark Street  
Suite 1425  
Chicago IL 60654  
T +1 312 585 8329  
F +1 312 644 0738  
info@hawo-usa.com  
www.hawo-usa.com

**hawo ASIA**  
25 International Business Park  
#03-105 German Centre  
Singapore 609916  
T +65 6433 5339  
F +65 6433 5359  
info@hawo-asia.com  
www.hawo-asia.com





***LABORATORY EQUIPMENT MAINTENANCE, REPAIR, CALIBRATION AND SALES***

Established in 1987, Henderson Biomedical is the UK's leading laboratory equipment sales and service provider. Our knowledgeable team can provide you with excellent sales advice on a range of different types of laboratory equipment including centrifuges, refrigerators, freezers and heat sealers.

Henderson Biomedical is also able to provide you with first class after-sales service and calibration of your laboratory equipment. We are an **ISO 17025 (UKAS) accredited calibration test laboratory** and our team of Field Service Engineers cover the whole of the United Kingdom.

Please contact us for more information on the types of equipment we supply and the different after-sales services we can offer.

**Henderson Biomedical**

Unit 3, Swan Close  
Croydon CR0 2DZ  
United Kingdom

**Tel: 020 8663 4610**

**For sales enquiries: [sales@henderson-biomedical.co.uk](mailto:sales@henderson-biomedical.co.uk)  
For all other enquiries: [info@henderson-biomedical.co.uk](mailto:info@henderson-biomedical.co.uk)  
[www.henderson-biomedical.co.uk](http://www.henderson-biomedical.co.uk)**

