EBA 280 / 280 S

Inhalt des Dokuments / content of the document

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Operating instructions (EN)

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Operating instructions

EBA 280 / 280 S

Translation of the original operating instructions

©2023 – All rights reserved

Andreas Hettich GmbH & Co. KG

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D-78532 Tuttlingen, Germany

Telephone: +49 (0)7461 705-0

Fax: +49 (0)7461 705-1125

Email: info@hettichlab.com, service@hettichlab.com

Internet: www.hettichlab.com

Table of contents

1	About this document.	5
1.1	Use of this document.	5
1.2	Gender reference.	5
1.3	Symbols and labels in this document.	5
2	Safety.	5
2.1	Intended use.	5
2.2	Personnel requirements.	6
2.3	Operator's responsibility.	7
2.4	Safety instructions.	7
3	Device overview.	9
3.1	Technical data.	9
3.2	European registration.	11
3.3	Important labels on the packaging.	11
3.4	Important labels on the device.	12
3.5	Operating and indicator elements.	13
3.5.1	Control.	13
3.5.2	Indicator elements.	13
3.5.3	Controls.	13
3.6	Original spare parts.	14
3.7	Scope of delivery.	14
3.8	Returns.	15
4	Transport and storage.	15
4.1	Transport and storage conditions.	15
4.2	Fastening the transport lock.	16
5	Commissioning.	16
5.1	Unpacking the centrifuge.	16
5.2	Removing the transport lock.	17
5.3	Setting up and connecting the centrifuge.	17
5.4	Switching the centrifuge on and off.	19
6	Operation	19
6.1	Opening and closing the lid.	19
6.2	Removing and installing the rotor.	20
6.3	Loading.	21
6.4	Centrifugation.	23
6.4.1	Centrifugation in continuous operation.	23
6.4.2	Centrifugation with time preselection.	23
6.4.3	Short-time centrifugation.	24
6.4.4	Changing settings during centrifugation.	24
6.5	Quick stop function.	24
7	Software operation.	24
7.1	Centrifugation parameters.	24
7.1.1	Ramp-up and ramp-down parameters.	24
7.1.2	Runtime.	25

7.1.3	Speed, RPM.	25
7.1.4	Relative centrifugal force, RCF.	26
7.1.5	Relative centrifugal force RCF and centrifuging radius RAD	26
7.1.6	Centrifugation of substances or mixtures of substances with a density higher than 1.2 kg/dm ³	27
7.2	Programming.	27
7.2.1	Opening or loading programs.	27
7.2.2	Entering or changing programs.	28
7.2.3	Saving programs.	28
7.2.4	Automatic buffer.	28
7.3	Rotor detection.	29
7.4	Machine Menu.	29
7.4.1	Querying system information.	29
7.4.2	Cycle counter.	30
7.4.3	Querying operating hours and centrifugation runs.	31
7.4.4	Enabling or disabling dual time.	31
7.4.5	Audible signal.	32
7.4.5.1	General.	32
7.4.5.2	Setting an audible signal.	32
7.4.6	Visual signal.	33
7.4.7	Automatic unlocking of the lid.	33
7.4.8	Centrifugation data displayed after switching on.	34
7.4.9	Indicator backlight.	34
8	Cleaning and care.	35
8.1	Overview table.	35
8.2	Instructions for cleaning and disinfection.	35
8.3	Cleaning.	36
8.4	Disinfection.	36
8.5	Maintenance.	37
9	Troubleshooting.	38
9.1	Fault description.	38
9.2	Performing a MAINS RESET.	40
9.3	Emergency release.	40
9.4	Replacing the mains input fuse.	40
10	Disposal.	41
10.1	General instructions.	41
11	Index.	43

1 About this document

1.1 Use of this document

- Read this document carefully and in full before commencing initial operation of the device.
Observe other enclosed instruction sheets where necessary.
- This document constitutes part of the device and must be kept within easy reach.
- This document must be included if the device is passed on to a third party.
- The most recently updated version of this document in the available languages can be found on the manufacturer's website at: ➔ <https://www.hettichlab.com/de/download-center/>

1.2 Gender reference

The employed masculine or feminine language form is to facilitate reading. In the spirit of equal treatment, corresponding terms apply in principle to all genders and do not imply any valuation.

1.3 Symbols and labels in this document

General symbols

The following markers are used in this document to highlight instructions, results, listings, references and other elements:

Marker	Explanation
1. ➔ 2. ➔ 3. ➔ ... ➔	Step-by-step instructions
➔	Results of action steps
➔	References to sections of the document and other applicable documents
■ ... ■ ...	Listings without a fixed order
[Buttons]	Controls (for example: buttons, switches)
'Indicator'	Indicator elements (for example: signal lights, screen elements)

2 Safety

2.1 Intended use

Intended use

The centrifuge **EBA 280 / 280 S** is an in vitro diagnostic medical device in accordance with the In Vitro Diagnostic Medical Devices Regulation (EU) 2017/746. The device is used for centrifugation as well as enrichment of sample material of human origin for subsequent further processing for diagnostic purposes. The user can set each of the variable physical parameters within the limits set by the device.

The centrifuge may only be used by qualified personnel in closed laboratories. The centrifuge is only intended for the use referred to above. Intended use also includes observing all instructions in the user manual and compliance with inspection and maintenance. Any other use or use beyond this is considered improper. Andreas Hettich GmbH & Co. KG shall not be liable for any damage arising from this.

Non-intended use

- The centrifuge is not suitable for use in explosive or radioactive, or biologically or chemically-contaminated atmospheres.
- The user must take appropriate actions when centrifuging hazardous substances or mixtures of substances that are toxic, radioactive or contaminated with pathogenic microorganisms.
The manufacturer generally recommends using only centrifuge tubes with special screw caps designed for use with hazardous substances. Use sealable centrifuge tubes with a biosafety system for materials of risk groups 3 and 4.
- The manufacturer does not recommend centrifugation of flammable or explosive materials.
- The manufacturer does not recommend centrifugation of materials that react chemically with one another with high activation energy.

Foreseeable misuse

The manufacturer recommends using only accessories approved that it has approved for the intended purpose.

Only operate the centrifuge under supervision.

2.2 Personnel requirements

Required qualifications

The user has read the Operating Manual in full and familiarised themselves with the device.

NOTICE

Damage to the device by unauthorised personnel

- Tampering with and modifications to devices by unauthorised persons are at the operating organisation's own risk and will result in the loss of all warranty and liability claims.

Trained user

The user has been educated and trained in laboratory work and is able to carry out the work assigned to them, and to recognise and prevent potential hazards independently.

Personal protective equipment

Lack of personal protective equipment or unsuitable personal protective equipment increases the risk of impaired health and injury.

- Only use personal protective equipment that is in proper condition.
- Only use personal protective equipment that is adapted to the person (correct size, for example).
- Observe instructions on other protective equipment for specific activities.

2.3 Operator's responsibility

Follow the instructions in this document for proper and safe use of the device.

Keep the user manual for future reference.

Provide information

- Following the instructions in this document will help:
 - To avoid dangerous situations.
 - To minimise repair costs and downtime.
 - To increase the reliability and service life of the device.
- The operator is responsible for compliance with company regulations, standards and national laws.
- Note and keep the revision of the document separate from the document. If lost, the document can be replaced in the correct revision.
- Keep the user manual available at the place where the device is used.
- Pass the user manual on to the buyer when the device is sold.

Personnel training

Lack of knowledge when working with the device may result in serious injury or death.

- Instruct personnel on their tasks and the associated risks in accordance with the instruction.

2.4 Safety instructions

Reporting serious incidents and notifiable incidents

In the event of serious incidents or notifiable incidents involving the device or its accessories, these must be reported to the manufacturer and, where applicable, to the competent authority where the user and/or the patient is registered.

DANGER

Risk of contamination for the user due to inadequate cleaning or failure to observe the cleaning instructions.

- Observe cleaning instructions.
- Wear personal protective equipment when cleaning the device.
- Observe laboratory regulations (e.g. TRBAs, the German Protection against Infection Act, hygiene plan) for handling biological agents.

DANGER

Fire and explosion hazard due to hazardous substances in samples.

- Observe relevant regulations and directives for handling chemicals and hazardous substances.
- Do not use aggressive chemicals (for example: dangerous, corrosive extraction agents such as chloroform, strong acids).

**WARNING**

Dangers due to insufficient maintenance or maintenance not carried out on time.

- Follow maintenance intervals.
- Check the device for visible damage or defects.
If any visible damage or defects are present, take the device out of service and inform a service technician.

**! WARNING**

Risk of electric shock due to ingress of water or other liquids.

- Protect the device against external liquids.
- Do not pour any liquids into the interior of the device.
- Transport using original transport packaging.

**! WARNING**

Contamination with hazardous substances and substance mixtures!

Observe the following actions for substances and substance mixtures that are toxic, radioactive and/or contaminated with pathogenic microorganisms:

- As a rule, use only centrifuge tubes with special screw caps for hazardous substances.
- Use sealable centrifuge tubes with a biosafety system for materials of risk groups 3 and 4.
- If no biosafety system is used, the device is not micro-biologically tight in the sense of standard EN / IEC 61010-2-020.
- Contact the manufacturer if necessary.

**CAUTION**

Risk of injury due to rotating rotor

Long hair and items of clothing can get caught on the rotor if the rotor is moved manually.

- Tie long hair back.
- Do not allow garments to hang in the centrifuging chamber.

**NOTICE**

Damage to the device electronics due to incorrect voltage or frequency at the device circuit breaker.

- Operate the device with the correct mains voltage and mains frequency.

The value can be found in the technical data and on the rating plate.

NOTICE

Damage to the device and samples due to premature program termination.

Premature program termination is caused by power failure, switching off during the program or pulling out the mains plug.

- Do not switch off the device while the program is running.
- Do not trigger the emergency release on the device while the program is running.
- Do not pull out the mains plug while the program is running.

3 Device overview

3.1 Technical data

Manufacturer	Andreas Hettich GmbH & Co. KG, D-78532 Tuttlingen			
Model	EBA 280		EBA 280 S	
Type	1101	1101-01	1102	1102-01
Mains voltage ($\pm 10\%$)	200-240 V 1~	100-127 V 1~	200-240 V 1~	100-127 V 1~
Mains frequency	50-60 Hz	50-60 Hz	50-60 Hz	50-60 Hz
power consumption	185 VA	185 VA	330 VA	330 VA
Power consumption	0.85 A	1.75 A	1.6 A	3.0 A
max. capacity	6 x 50 ml			
max. permissible density	1.2 kg/dm ³			
max. speed (RPM)	6000			
max. acceleration (RCF)	4146		5071	
max. kinetic energy	2700 Nm			
Obligation to perform checks (DGUV Rules 100-500) (valid only in Germany)	No			
Ambient conditions (EN / IEC 61010-1):				
Installation site	indoors only			
Altitude	up to 2000 m above sea level			
Ambient temperature	2 °C to 40 °C			

Humidity	maximum relative humidity 80% for temperatures up to 31 °C, decreasing linearly to 50% relative humidity at 40 °C.		
Overvoltage category (IEC 60364-4-443)	II		
Pollution level	2		
Device protection class	I not suitable for use in potentially explosive atmospheres.		
EMC:			
Emitted EM interference, EM interference immunity	EN / IEC 61326-1 Class B	FCC Class B	EN / IEC 61326-1 Class B FCC Class B
Noise level (rotor-dependent)	≤51 dB(A)		≤56 dB(A)
Dimensions:			
Width	326 mm		
Depth	400 mm		
Altitude	242 mm		
Weight	approx. 11 kg		

Rating plate


Fig. 1: Rating plate

- 1 Item number
- 2 Serial number
- 3 Revision
- 4 Equipment number
- 5 Data matrix code
- 6 any labelling indicating whether medical device or in vitro diagnostic medical device
- 7 Global Trade Item Number (GTIN)
- 8 Date of manufacture
- 9 Serial number
- 10 any EAC mark, CE mark

- 11 Country of manufacture
- 12 Date of manufacture
- 13 Mains frequency
- 14 Maximum kinetic energy
- 15 Maximum permissible density
- 16 Manufacturer's address
- 17 any Coolant circuit pressure
- 18 any Coolant capacity
- 19 any Coolant type
- 20 Revs per minute
- 21 Performance values
- 22 Mains voltage
- 23 any Device designation
- 24 Manufacturer's logo

3.2 European registration

Device conformity

Device conformity according to EU directives.

Single Registration Number

SRN: DE-MF-000010680

Basic-UDI-DI

Basic-UDI-DI	Device assignment
040506740100089Y	EBA 280 / 280 S (in vitro diagnostic medical device)

3.3 Important labels on the packaging

TOP

This is the correct upright position of the shipping container for transport and/or storage.

FRAGILE GOODS

The contents of the shipping container are fragile, so it must be handled with care.

PROTECT FROM MOISTURE

The shipping container must not be exposed to rain and kept in a dry environment.

TEMPERATURE LIMITATION
The shipping container must be stored, transported and handled within the indicated temperature range (-20 °C to +60 °C).

HUMIDITY LIMITATION
The shipping container must be stored, transported and handled within the indicated humidity range (10% to 80%).

STACK LIMITATION BASED ON QUANTITY
Maximum number of identical packages that may be stacked on the lowest package, "n" standing for the number of packages allowed. The lowest package is not included in "n".

3.4 Important labels on the device

The signs on the device must not be removed or covered, or have anything pasted over them.

Attention, general danger area.
Ensure you read the instructions for commissioning and operation and observe the safety instructions before using the device.

Biohazard warning.

Direction of rotation of the rotor.
The orientation of the arrow indicates the rotor's direction of rotation.

Direction of rotation of the emergency release.

Symbol for the separate collection of electrical and electronic equipment, in accordance with Directive 2012/19/EU (WEEE).
Use in European Union countries, Norway and Switzerland.

3.5 Operating and indicator elements

3.5.1 Control


Fig. 2: Control

3.5.2 Indicator elements

Fig. 3: 'Lid unlocked' indicator

- The indicator appears when the lid is unlocked.

Fig. 4: 'Lid locked' indicator

- The indicator appears when the lid is locked.

Fig. 5: 'Rotation' indicator

- The indicator light rotates when the rotor is turning.

Fig. 6: [RCF] button

- The indicator above the button lights up when the RCF is displayed.

Fig. 7: [START/PULSE] button

- The button lights up during the centrifugation run for as long as the rotor is not yet at a standstill.

Fig. 8: [STOP/OPEN] button

- The right side of the button lights up when the centrifuge is in ramp-down. The rotor has not yet stopped.
- The left side of the button lights up when the rotor is stationary.
- The light on the left side of the button goes out when the lid is unlocked.

3.5.3 Controls

Fig. 9: [Mains switch]

- Switch the device on and off.

Fig. 10: [PROG] button

- Retrieve programs.
- Save programs.


Fig. 11: [RCF] button

- Toggle between RCF indicator and RPM indicator.
- Relative centrifugal force, RCF. The RCF is displayed in brackets > <.
- Speed, RPM.

Fig. 12: [SELECT] button

- Selecting the individual parameters.
- Open 'MACHINE MENU'.
- Scroll forward in the menus.

Fig. 13: [START/PULSE] button

- Start centrifugation run.
- Short-time centrifugation. The centrifugation run takes place as long as the button is being pressed.
- Save entries and changes.

Fig. 14: [STOP/OPEN] button

- End the centrifugation run. The rotor coasts to a stop with the preselected ramp-down parameter.
- Pressing the button twice triggers the quick stop function.
- Unlock the lid.
- Exit parameter input and the menus

Fig. 15: Configuration buttons

- Change the value of a parameter.
- The value changes at an increasing rate if the button is held down.

3.6 Original spare parts

Use only original spare parts from the manufacturer and approved accessories.

3.7 Scope of delivery

The following accessories are supplied with the centrifuge:

- 2 Fuse link
- 1 hex key (SW5 x 100)
- 1 power cable
- 1 user manual
- 1 instruction sheet, transport lock
- 1 Emergency release instruction sheet

Rotors and the corresponding accessories are supplied depending on the order.

3.8 Returns

An original Return Material Authorisation (RMA) form from the manufacturer must always be requested for a return. Secure and reliable acceptance and booking in of the goods with the manufacturer is not possible without an original RMA form from the manufacturer. The Return Material Authorisation (RMA) form contains a Declaration of No Objection (UBE), which must be completed in full and enclosed with the return.

If the device and/or accessories are returned to the manufacturer, the complete return shipment must be cleaned and decontaminated by the sender. If returns are not cleaned and/or decontaminated or are insufficiently cleaned and/or decontaminated, this will be performed by the manufacturer and charged to the sender.

The original transport locks must be attached for return shipment, see [Chapter 4 'Transport and storage' on page 15](#). The device must be shipped in its original packaging.

4 Transport and storage

4.1 Transport and storage conditions

Transport conditions

NOTICE

Damage to the device due to failure to use the transport locks.

- Secure the transport locks before transporting the device.

NOTICE

Damage to the device due to condensation.

There is a risk of condensation forming on electrical components when component surfaces are cold and the surrounding air is warmer. The condensation that forms may cause a short circuit and/or destroy electronics.

- Warm the device up for at least 3 hours in a warm room before connecting it to the mains.
or
- Warm up for 30 minutes in a cold room.

- Before transporting, fasten the transport lock and disconnect the device from the mains socket.
- The transport temperature must be between -20 °C and +60 °C.
- Humidity must not be condensing. Humidity must be between 10% and 80%.
- Be aware of the weight of the device.
- When transporting using a transport aid (e.g., a pallet truck), the transport aid must be able to carry at least 1.6 times the transport weight of the device.
- Secure the device to prevent it tipping over and falling down during transport.
- Never transport the device sideways or upside down.

Storage conditions

- The device must be stored in the original packaging.
- Only store the device in dry rooms.

- The storage temperature must be between -20 °C and +60 °C.
- Humidity must not be condensing. Humidity must be between 10% and 80%.

4.2 Fastening the transport lock

Personnel:

- Trained user

The lid is closed.

The mains cable is disconnected from the device.

1. ➤ Tilt the device on the back of the device.
2. ➤ Insert 2 spacer sleeves (1).
3. ➤ Screw in 2 screws (2).


Fig. 16: Transport lock

- 1 Spacer sleeves
2 Screws

5 Commissioning

5.1 Unpacking the centrifuge

CAUTION

Danger of crushing due to parts falling out of the transport packaging.

- Keep the device balanced during the unpacking process.
- Only open the packaging at the points provided for this purpose.

CAUTION

Risk of injury from lifting heavy loads.

- Provide an adequate number of helpers.
- Note the weight. See ➔ Chapter 3.1 'Technical data' on page 9.

NOTICE

Damage to the device due to improper lifting.

- Do not lift the centrifuge by the control panel or the control panel holder.

Personnel:

- Trained user

1. ➤ Open the box at the top.

2. ➤ Remove the padding.
3. ➤ Remove the device and accessories by lifting them up out of the box.
4. ➤ Place the device on a stable and level surface.

5.2 Removing the transport lock

Personnel:

- Trained user

The lid is closed.

The mains cable is disconnected from the device.

1. ➤ Tilt the device on the back of the device.
2. ➤ Unscrew 2 screws (2).
3. ➤ Remove 2 spacer sleeves (1).
4. ➤ Keep the screws and spacer sleeves in a safe place.


Fig. 17: Transport lock

- 1 Spacer sleeve
- 2 Screw

5.3 Setting up and connecting the centrifuge

Setting up the centrifuge

WARNING

Risk of injury due to failing to maintain a sufficient distance to the centrifuge.

- As per EN / IEC 61010-2-020, no persons, hazardous materials or objects may be present within a **safety zone of 300 mm** around the centrifuge during a centrifugation run.
- A distance of **300 mm** from the ventilation slots and ventilation openings of the centrifuge must be maintained.

CAUTION

Risk of crushing and damage to the device due to it falling down because of vibration-induced position alterations.

- Place the device on a stable and level surface.
- Select the installation surface dependent on the weight of the device.

**NOTICE**

Damage to the samples and the device if the ambient temperature exceeds or falls below the respective maximum/minimum permissible ambient temperature.

- Comply with the maximum and minimum permissible ambient temperatures for installation of the device.
- Do not place the device next to a heat source.
- Do not expose the device to direct sunlight.
- Do not expose the device to frost.

Personnel:

- Trained user

1. Place the device on a stable and level surface.
2. Maintain a distance of 300 mm around the device.
3. Comply with the ambient conditions in the technical data (→ *Chapter 3.1 'Technical data' on page 9*).

Connecting the centrifuge**NOTICE**

Damage to the device by unauthorised personnel

- Tampering with and modifications to devices by unauthorised persons are at the operating organisation's own risk and will result in the loss of all warranty and liability claims.

**NOTICE**

Damage to the device due to condensation.

There is a risk of condensation forming on electrical components when component surfaces are cold and the surrounding air is warmer. The condensation that forms may cause a short circuit and/or destroy electronics.

- Warm the device up for at least 3 hours in a warm room before connecting it to the mains.
- or
- Warm up for 30 minutes in a cold room.

Personnel:

- Trained user

1. A type B residual current circuit breaker must be used if the device is additionally protected with a residual current circuit breaker in the building installation.

When using a different type, the residual current circuit breaker may either not switch off the unit if there is a fault on the unit, or it may switch off the unit even though there is no fault on the unit.

2. Check whether the mains voltage matches the specification on the rating plate.
3. Connect the device to a standard mains socket using the mains cable.

5.4 Switching the centrifuge on and off.

Switching the centrifuge on

Personnel:

- Trained user

→ Set the mains switch to [I].

- ➔ The buttons flash, depending on the centrifuge type.

The following indicators appear one after the other, depending on the centrifuge type:

- the centrifuge model
- the machine type and program version
- the rotor code (R) and the maximum speed (maxRPM) of the last rotor recognised by the rotor detection
- When the lid is closed: For centrifuges with cooling, 'Open the lid' indicator. The lid opens in centrifuges without cooling.
- When the lid is open: The centrifugation data of the last program used or program 1.

Switching off the centrifuge

The rotor is stationary.

→ Set the mains switch to [0].

6 Operation

6.1 Opening and closing the lid

Opening the lid

Personnel:

- Trained user

The centrifuge is switched on.

The rotor is stationary.

→ Press the [STOP/OPEN] button.

- ➔ The lid unlocks by means of a motor.

The 'Lid unlocked' indicator appears.

Closing the lid

⚠ CAUTION

Crushing hazard when closing the lid.

Danger of fingers getting crushed when the closing motor pulls the lid against the seal.

- No parts of the operator's body should be in the hazard zone of the lid when closing the lid.
- To close the lid, press on the lid from above.

NOTICE

Damage to the device caused by the lid slamming.

- Close the lid slowly.
- Do not slam the lid.

Personnel:

- Trained user

- Close the lid and press the front edge of the lid down gently.
 - ➔ The lid locks using a motor.
 - The 'Lid locked' indicator appears.

6.2 Removing and installing the rotor

Removing a rotor with a release knob


Fig. 18: Rotor installation and removal

- 1 Release knob
- 2 Motor shaft

Personnel:

- Trained user

1. → Open the lid.
2. → Pull the release knob (1) up, hold it and lift the rotor off the motor shaft (2).

Installing a rotor with a release knob

Personnel:

- Trained user

The lid is open.

1. → Clean the motor shaft (2) and rotor hole.
2. → Lightly grease the motor shaft (2), see → Chapter 8.2 'Instructions for cleaning and disinfection' on page 35.
3. → Lift the rotor by the release knob (1) and place it vertically on the motor shaft (2).
 - ➔ The rotor engages automatically on the motor shaft.
4. → Check that the rotor is firmly seated by holding the rotor on the left and right and pulling up gently.

Removing a rotor without a release knob


Fig. 19: Rotor installation and removal

- 1 Handle
- 2 Motor shaft

Personnel:

- Trained user

1. → Open the lid.
2. → Hold the rotor by the handle (1) and lift it off the motor shaft (2).

Installing a rotor without a release knob

Personnel:

- Trained user
1. → Clean the motor shaft (2) and rotor hole.
 2. → Lightly grease the motor shaft (2).
 3. → Hold the rotor by the handle (1), place it vertically on the motor shaft (2) and press it down as far as it will go.

6.3 Loading

Filling centrifuge tubes

WARNING

Risk of injury from contaminated sample material.

Contaminated sample material escapes from the sample tube during centrifugation.

- Use centrifuge tubes with special screw caps for hazardous substances.
- For risk group 3 and 4 materials, use a biosafety system in addition to the sealable centrifuge tubes (see WHO's 'Laboratory Biosafety Manual').

NOTICE

Damage to the device due to highly corrosive substances.

Highly corrosive substances may impair the mechanical strength of rotors, buckets and accessories.

- Do not centrifuge highly corrosive substances.

Standard glass centrifuge tubes can be loaded up to RCF 4000 (DIN 58970 part 2).

Personnel:

- Trained user
- Fill centrifuge tubes outside the centrifuge.

The maximum capacity of the centrifuge tubes specified by the manufacturer must not be exceeded.





With angle rotors, the centrifuge tubes must only be filled to the extent that no liquid can be ejected from the tubes during the centrifugation run.

It must be ensured that there is a uniform fill level in the tubes in order to keep the weight differences in the centrifuge tubes as low as possible.

Loading swing-out rotors

Personnel:

- Trained user

1. → Check that the rotor is firmly seated.
2. → The centrifuge tubes must be distributed symmetrically across all rotor locations.

The weight of the permissible filling capacity is indicated on each rotor. The weight must not be exceeded.

No liquid must be allowed to enter the buckets and the centrifuging chamber when loading the buckets and swinging them out during the centrifugation run.

For containers with rubber inserts, there must always be the same number of rubber inserts under the centrifuge tubes.

All rotor locations must be filled with the same buckets. Certain buckets are marked with the number of the rotor location. The buckets must only be inserted in the corresponding rotor location.

Buckets marked with a set number (for example S001/4) must only be used in the set.

Loading the angle rotors

Personnel:





- Trained user

1. → Check that the rotor is firmly seated.
2. → The centrifuge tubes must be distributed evenly over all locations on the rotor.

No liquid must be allowed to enter the rotor and the centrifuging chamber when loading the rotor.

With rotors, the centrifuge tubes must only be filled to the extent that no liquid can be ejected from the tubes during the centrifugation run.

The weight of the permissible filling capacity is indicated on each rotor. The weight must not be exceeded.

6.4 Centrifugation

6.4.1 Centrifugation in continuous operation

Personnel:

- Trained user

1.  If required: Press the *[RCF]* button.
 - The parameter RCF (*'>RCF<'*) or RPM (*'RPM'*) is displayed. Press the *[RCF]* button to toggle between the two parameters.
2.  Enter the desired speed (RPM) or relative centrifugal force (RCF).
See ➔ *Chapter 7.2.1 'Opening or loading programs' on page 27.*
3.  Set the parameters t/min and t/sec to zero.
 - *'--:--'* is displayed.
See ➔ *Chapter 7.2.1 'Opening or loading programs' on page 27.*
4.  Press the *[START/PULSE]* button.
 - The centrifugation run is started.
The timing starts at *'0:00'*.
The rotor speed or the resulting RCF value and the elapsed time are displayed during the centrifugation run.
5.  Press the *[STOP/OPEN]* button to cancel the centrifugation run.
 - Ramp-down takes place with the set brake level. The brake level is displayed.
When the rotor is at a standstill, the lid opens, an audible signal sounds and the remaining number of run cycles (centrifugation runs) is displayed.

6.4.2 Centrifugation with time preselection

Personnel:



- Trained user

1.  If required: Press the *[RCF]* button.
 - The parameter RCF (*'>RCF<'*) or RPM (*'RPM'*) is displayed. Press the *[RCF]* button to toggle between the two parameters.
2.  Enter the desired speed (RPM) or relative centrifugal force (RCF).
See ➔ *Chapter 7.2.1 'Opening or loading programs' on page 27.*
3.  Set the parameters t/min and t/sec to the desired value.
See ➔ *Chapter 7.2.1 'Opening or loading programs' on page 27.*
4.  Press the *[START/PULSE]* button.
 - The centrifugation run is started.
The rotor speed or the resulting RCF value and the remaining time are displayed during the centrifugation run.
5.  Press the *[STOP/OPEN]* button to cancel the centrifugation run.
or
Wait for the centrifugation time to elapse.
 - Ramp-down takes place with the set brake level. The brake level is displayed.
When the rotor is at a standstill, the lid opens, an audible signal sounds and the remaining number of run cycles (centrifugation runs) is displayed.

6.4.3 Short-time centrifugation

Personnel:

- Trained user

1.  If required: Press the *[RCF]* button.
 - ➔ The parameter RCF ('>RCF<') or RPM ('RPM') is displayed. Press the *[RCF]* button to toggle between the two parameters.
2.  Enter the desired centrifugation parameters.
See → Chapter 7.2.1 'Opening or loading programs' on page 27.
3.  Press and hold the *[START/PULSE]* button.
 - ➔ The centrifugation run is started.
The timing starts at '0:00'.
The rotor speed or the resulting RCF value and the elapsed time are displayed during the centrifugation run.
4.  Release the *[START/PULSE]* button to end the centrifugation run.
 - ➔ Ramp-down takes place with the set brake level. The brake level is displayed.
When the rotor is at a standstill, the lid opens, an audible signal sounds and the remaining number of run cycles (centrifugation runs) is displayed.

6.4.4 Changing settings during centrifugation

The runtime, speed, relative centrifugal force (RCF) and ramp-up and ramp-down parameters can be changed during centrifugation.

-  Change the value of the desired parameter.
 - ➔ The values of the current program are copied to program location '#' and updated with the changed value.
The original program is not overwritten.

6.5 Quick stop function

Personnel:

- Trained user


-  Press the *[STOP/OPEN]* button twice.
 - ➔ Ramp-down with brake level "9" (shortest ramp-down time) is displayed and executed.

7 Software operation

7.1 Centrifugation parameters

7.1.1 Ramp-up and ramp-down parameters

Ramp-up level

1.  Press the *[SELECT]* button repeatedly until  is displayed.
2.  Use the *[Adjustment buttons]* to set the desired value.
A numerical value can be set from 1 to 9.
Adjustable in 1 second increments.
9 = shortest ramp-up time
1 = longest ramp-up time

3. Press the *[START/PULSE]* button.
 - The setting is transferred to the indicator.

Brake level

1. Press the *[SELECT]* button repeatedly until is displayed.
2. Use the *[Adjustment buttons]* to set the desired value.
A numerical value can be set from 0 to 9.
Adjustable in 1 second increments.
9 = shortest ramp-down time
1 = long ramp-down time
0 = longest ramp-down time (unbraked ramp-down).
3. Press the *[START/PULSE]* button.
 - The setting is transferred to the indicator.

7.1.2 Runtime Changing the runtime

For continuous operation, the minutes and seconds must be set to zero.

1. Press the *[SELECT]* button repeatedly until *'t/min'* is displayed.
2. Use the *[Adjustment buttons]* to set the desired value.
A numerical value can be set from 1 to 99 minutes.
Adjustable in 1 minute increments.
3. Press the *[SELECT]* button.
 - *'t/sec'* is displayed.
4. Use the *[Adjustment buttons]* to set the desired value.
A numerical value can be set from 1 to 59 seconds.
Adjustable in 1 second increments.
5. Press the *[START/PULSE]* button.
 - The settings are shown in the indicator.

Start of runtime count

- The *'Dual time'* function is enabled. The function is enabled ex works.
1. Press the *[SELECT]* button repeatedly until *'Begins at START'* or *'Begins at SPEED'* is displayed.
 2. Use the *[Adjustment buttons]* to select the desired setting.
 - *'Begins at START'* = Runtime counting begins after the start of the centrifugation run.
 - *'Begins at SPEED'* = Timing of the runtime starts after the set speed is reached.
This is indicated by the symbol in the indicator to the left of the time.
 3. Press the *[START/PULSE]* button.
 - The settings are shown in the indicator.

7.1.3 Speed, RPM

1. Press the *[SELECT]* button repeatedly until *'RPM'* is displayed.

2. Use the *[Adjustment buttons]* to set the desired value.
A numerical value from 200 RPM to the maximum rotor speed can be set.
Adjustable in 10 second increments.
3. Press the *[START/PULSE]* button.
➔ The setting is transferred to the indicator.

7.1.4 Relative centrifugal force, RCF

The relative centrifugal force RCF is dependent on the speed and the centrifuging radius.

The relative centrifugal force RCF is stated as a multiple of the acceleration due to gravity (g).

The relative centrifugal force RCF is a dimensionless numerical value and is used to compare the separation and sedimentation performance.

$$RCF = \left(\frac{RPM}{1000} \right)^2 * r * 1,118$$

$$RPM = \sqrt{\frac{RCF}{r * 1,118} * 1000}$$

RCF = Relative Centrifugal Force

RPM = speed

r = centrifuging radius in mm = distance from the centre of the axis of rotation to the bottom of the centrifuge tube.

7.1.5 Relative centrifugal force RCF and centrifuging radius RAD

The relative centrifugal force RCF is dependent on the centrifuging radius RAD. The centrifuging radius must be set before setting the relative centrifugal force.

1. Press the *[RCF]* button.
➔ The LED above the button lights up.
2. Press the *[SELECT]* button repeatedly until 'R:' and 'RCF' are displayed.
➔ Value of the 'RCF' parameter is displayed in brackets < >.
3. Use the *[Adjustment buttons]* to set the desired 'RCF'.
A numerical value can be set that gives a speed between 200 RPM and the maximum rotor speed.
Adjustable in 1 second increments.
The set centrifuging radius is displayed while entering the RCF.
4. Press the *[SELECT]* button repeatedly until 'RAD/mm' is displayed.
5. Use the *[Adjustment buttons]* to set the desired centrifuging radius.
A numerical value from 10 mm to 330 mm can be set.
Adjustable in 1 millimetre increments
The value of the RCF adjusts automatically when the centrifuging radius is changed.
6. Press the *[START/PULSE]* button.
➔ The setting is transferred to the indicator.

7.1.6 Centrifugation of substances or mixtures of substances with a density higher than 1.2 kg/dm³

The density of the substances or mixtures of substances must not exceed 1.2 kg/dm³ during centrifugation at maximum speed. The speed must be reduced for substances or substance mixtures with a higher density. The permissible speed can be calculated using the following formula:

$$\text{Reduced speed } (n_{red}) = \sqrt{\frac{1,2}{\text{Greater density [kg/dm}^3\text{]}} * \text{maximum speed [RPM]}}$$

For example: Maximum speed 4000 RPM, density 1.6 kg/dm³

$$n_{red} = \sqrt{\frac{1,2(\text{kg/dm}^3)}{1,6(\text{kg/dm}^3)}} * 4000 \text{ RPM} = 3464 \text{ RPM}$$

If, in exceptional cases, the maximum load indicated on the bucket is exceeded, the speed must also be reduced. The permissible speed can be calculated using the following formula:

$$\text{Reduced speed } (n_{red}) = \sqrt{\frac{\text{maximum load [g]}}{\text{actual load [g]}}} * \text{maximum speed [RPM]}$$

For example: Maximum speed 4000 RPM, maximum load 300 g, actual load 350 g

$$n_{red} = \sqrt{\frac{300 \text{ g}}{350 \text{ g}}} * 4000 \text{ RPM} = 3703 \text{ RPM}$$

Please contact the manufacturer if you are not sure.

7.2 Programming

7.2.1 Opening or loading programs

1.  Use the *[PROG]* button to select the 'PROG RCL' parameter.
2.  Use the *[Adjustment buttons]* to set the desired program location.
Program locations 1 to 9 and # can be set.
3.  Press the *[START/PULSE]* button.
 - ◆ '*** OK ***' is displayed briefly.
 - The centrifugation data of the desired program location is displayed
4.  To check the parameters: Press the *[SELECT]* button.
5.  To exit the parameter indicator: Press the *[OPEN/STOP]* button or do not press any button for 8 seconds.

7.2.2 Entering or changing programs


Program locations 1-9 are available.

No programs may be stored in program location #. Program location # serves as a buffer for changed centrifugation parameters.


A dash "-" appears in the indicator instead of the program location number if centrifugation parameters are changed but then not saved to a program location. The centrifugation parameters are automatically saved in program location # after starting the centrifugation run.

The centrifugation parameters in program location # are overwritten every time a centrifugation run is executed with changed centrifugation parameters that have not been saved to a program location.

1.  If required: Press the *[RCF]* button to toggle between RPM and RCF indicator.
 - The indicator above the button lights up.
2.  If required: Press the *[SELECT]* button to select the desired parameter and set it with the *[Adjustment buttons]*.
The parameters t/min and t/sec must be set to 0 with the *[Adjustment buttons]* to set continuous operation. Continuous operation is shown in the time indicator with '--:--'.
3.  Use the *[SELECT]* button to select the 'PROG STO' parameter.
4.  Use the *[Adjustment buttons]* to set the desired program location.
5.  Press the *[START/PULSE]* button.
 - Settings are stored in the desired program location.
**** OK **** is displayed briefly.

The settings are always stored in program location # if the *[START/PULSE]* button is pressed without the 'PROG STO' parameter being selected.

7.2.3 Saving programs

1.  Press the *[PROG]* button twice.
 - 'PROG STO' is displayed.
PROG STO: Program location where the centrifugation parameters are saved.
2.  Use the *[Adjustment buttons]* to set the desired program location.
3.  Press the *[START/PULSE]* button.
 - Settings are stored in the desired program location.
**** OK **** is displayed briefly.



The settings are always stored in program location # if the *[START/PULSE]* button is pressed without the 'PROG STO' parameter being selected.

7.2.4 Automatic buffer

The centrifugation data is saved to program location '#' every time a centrifugation run is started.

No programs can be stored in program location '#'.

7.3 Rotor detection

- Rotor detection is performed after starting a centrifugation run.
 - The centrifugation run is cancelled after rotor detection if the rotor code of the rotor is being read in for the first time. The rotor code (R) and the newly detected rotor's maximum speed (maxRPM) permissible for the centrifuge are displayed.
1.  Press any button after the rotor has stopped.
 - The maximum permissible number of run cycles ('Cyc lim') is displayed.
 2.  Set the maximum permissible number of run cycles of the rotor or the bucket.
- Every rotor in a centrifuge has a maximum permissible speed. The maximum rotor speed ('ROTOR MAX') is displayed briefly after the centrifugation run has started. The centrifugation run is cancelled if the set speed is greater than the permissible maximum rotor speed. The permissible maximum rotor speed is displayed.
-  Adjust the speed to the maximum rotor speed.
- If the cycle counter is enabled, the remaining number of run cycles (centrifugation runs) of the rotor code used is displayed briefly after opening the lid.



7.4 Machine Menu

7.4.1 Querying system information

The following system information can be queried:

- Centrifuge model
- Maximum speeds of the various rotor codes
- Centrifuge program version
- Frequency converter type
- Program version for the frequency inverter

The rotor is stationary.

1.  Press and hold the [SELECT] button.
 - '*MACHINE MENU*' is displayed after 8 seconds.
2.  Press the [SELECT] button.
 - '-> Info' is displayed.
3.  Press the [START/PULSE] button.
 - The centrifuge model is displayed.
4.  Press the [SELECT] button.
 - The following is displayed:
 - 'R': Rotor code of the rotor
 - '*': The rotor code of the rotor currently being used is marked with an asterisk.
 - 'RPMmax': Maximum speed of the centrifuge rotor
5.  If necessary: press [Adjustment buttons].
 - The maximum speeds of the various rotor codes are displayed.
6.  Press the [SELECT] button.
 - The centrifuge program version 'CP FW=' is displayed.

7.  Press the *[SELECT]* button.
 - The type of frequency converter '*FC type LC*' of the centrifuge is displayed.
8.  Press the *[SELECT]* button.
 - The program version of the frequency converter '*FC FW=*' of the centrifuge is displayed.
9.  Press the *[STOP/OPEN]* button twice to exit the '*-> Info*' menu or
 Press the *[STOP/OPEN]* button three times to exit the '**MACHINE MENU**'.

7.4.2 Cycle counter

The centrifuge is equipped with a cycle counter. The cycle counter counts the run cycles (centrifugation cycles). The remaining number of run cycles (centrifugation runs) is displayed briefly after each centrifugation run.

The centrifugation run is cancelled when the rotor is recognised by the rotor detection for the first time. '*Cyc lim = (50000)*' is displayed after pressing any button. The maximum permissible number of run cycles indicated on the rotor must be entered before the centrifugation run can be restarted.

If the maximum permissible number of rotor run cycles entered is exceeded, '*Cycles passed*' is displayed after each start of a centrifugation run. The centrifugation run must be restarted. The rotor must be replaced with a new one.

Once the rotor has been replaced, the cycle counter must be reset to '0'.

Entering the maximum permissible number of run cycles


The maximum permissible number of run cycles must be entered after starting the first centrifugation run.

'*Cyc lim = (50000)*' is displayed.

1.  Use the *[Adjustment buttons]* to set the maximum permissible number of run cycles indicated on the rotor.
2.  Press the *[START/PULSE]* button.
 - The setting is stored.
 - '*Store cycles ...*' is displayed briefly.

Resetting the cycle counter

The cycle counter must be reset to '0' after installing a new rotor.

1.  Press and hold the *[SELECT]* button.
 - '**MACHINE MENU**' is displayed after 8 seconds.
2.  Press the *[SELECT]* button repeatedly until '*-> Time & Cycles*' is displayed.
3.  Press the *[START/PULSE]* button.
4.  Press the *[SELECT]* button repeatedly until '*Cyc sum=...*' is displayed.
5.  Press the *[RCF]* button.
6.  Press the *[▼]* button.
 - The number of run cycles completed is reset to zero.
7.  Press the *[START/PULSE]* button.
 - '*Store cycles...*' is displayed.

8.  Press the *[STOP/OPEN]* button twice to exit the '*-> Time & Cycles*' menu
or
Press the *[STOP/OPEN]* button three times to exit the *"*MACHINE MENU*"*.

7.4.3 Querying operating hours and centrifugation runs

The operating hours are divided into internal and external operating hours.

- Internal operating hours: Total time for which the device has been switched on.
- External operating hours: Total time of centrifugation runs to date.

The rotor is stationary.


1.  Press and hold the *[SELECT]* button.
 - ➔ *"*MACHINE MENU*"* is displayed after 8 seconds.
2.  Press the *[SELECT]* button repeatedly until '*-> Time & Cycles*' is displayed.
3.  Press the *[START/PULSE]* button.
 - ➔ '*TimeExt=*' is displayed.
TimeExt: External operating hours
4.  Press the *[SELECT]* button.
 - ➔ '*TimeInt=*' is displayed.
TimeInt: Internal operating hours
5.  Press the *[SELECT]* button.
 - ➔ '*Starts=*' is displayed.
Starts: Number of all centrifugation runs
6.  Press the *[STOP/OPEN]* button twice to exit the '*-> Time & Cycles*' menu
or
Press the *[STOP/OPEN]* button three times to exit the *"*MACHINE MENU*"*.

7.4.4 Enabling or disabling dual time

The '*Dual time*' function must be enabled to allow the '*Begins at SPEED*' parameter to be set. The function is enabled ex works.

The rotor is stationary.

1.  Press and hold the *[SELECT]* button.
 - ➔ *"*MACHINE MENU*"* is displayed after 8 seconds.
2.  Press the *[SELECT]* button repeatedly until '*-> Settings*' is displayed.
3.  Press the *[START/PULSE]* button.
 - ➔ '*End beep = on*' or '*End beep = off*' is displayed.
4.  Press the *[SELECT]* button repeatedly until '*Dual time = on*' or '*Dual time = off*' is displayed.
5.  Use *[Adjustment buttons]* to set '*off*' or '*on*'.
 - off = The function is disabled
 - on = The function is enabled.

6.  Press the *[START/PULSE]* button.
 - ➔ The settings are stored.
 - 'Store Settings...' is displayed briefly.
 - '-> Settings' is then displayed.
7.  Press the *[STOP/OPEN]* button once to exit the '-> Settings' menu or
Press the *[STOP/OPEN]* button twice to exit the '*MACHINE MENU*'.

7.4.5 Audible signal



7.4.5.1 General

The audible signal sounds:

- after a problem occurs in the 2 s interval.
- after completion of the centrifugation run and rotor standstill in the 30 s interval.

Opening the lid or pressing any button stops the audible signal.

7.4.5.2 Setting an audible signal

1.  Press and hold the *[SELECT]* button.
 - ➔ '*MACHINE MENU*' is displayed after 8 seconds.
2.  Press the *[SELECT]* button repeatedly until '-> Settings' is displayed.
3.  Press the *[START/PULSE]* button.
 - ➔ 'End beep = on' or 'End beep = off' is displayed.
4.  Use *[Adjustment buttons]* to set 'off' or 'on'.
 - off: Audible signal after completion of the centrifugation run is disabled.
 - on: Audible signal after completion of the centrifugation run is enabled.
5.  Press the *[SELECT]* button.
 - ➔ 'Error beep = on' or 'Error beep = off' is displayed.
6.  Use *[Adjustment buttons]* to set 'off' or 'on'.
 - off: Audible signal after the occurrence of a malfunction is disabled.
 - on: Audible signal after the occurrence of a malfunction is enabled.
7.  Press the *[SELECT]* button.
 - ➔ 'Beep volume = min', 'Beep volume = mid' or 'Beep volume = max' is displayed.
8.  Use *[Adjustment buttons]* to set 'min', 'mid' or 'max'.
 - min: The volume of the audible signal is set to low.
 - mid: The volume of the audible signal is set to medium.
 - Max: The volume of the audible signal is set to loud.
9.  Press the *[START/PULSE]* button.
 - ➔ The setting is stored.
 - 'Store Settings...' is displayed briefly.
 - '-> Settings' is then displayed.

10. Press the *[STOP/OPEN]* button once to exit the '-> Settings' menu
or
Press the *[STOP/OPEN]* button twice to exit the ***MACHINE MENU***.

7.4.6 Visual signal

The indicator backlight flashes as a visual signal after the centrifugation run is finished.

Switching on and off

1. Press and hold the *[SELECT]* button.
 - ***MACHINE MENU*** is displayed after 8 seconds.
2. Press the *[SELECT]* button repeatedly until '-> Settings' is displayed.
3. Press the *[START/PULSE]* button.
 - 'End beep = on' or 'End beep = off' is displayed.
4. Press the *[SELECT]* button repeatedly until 'End blinking=off' or 'End blinking =on' is displayed.
5. Use *[Adjustment buttons]* to set 'off' or 'on'.
 - off: Backlight does not flash.
 - on: Backlight flashes.
6. Press the *[START/PULSE]* button.
 - The setting is stored.
 - 'Store setting...' is displayed briefly.
 - '-> Settings' is then displayed.
7. Press the *[STOP/OPEN]* button once to exit the '-> Settings' menu
or
Press the *[STOP/OPEN]* button twice to exit the ***MACHINE MENU***.

7.4.7 Automatic unlocking of the lid

Setting whether or not the lid unlocks automatically after the centrifugation run.


The rotor is stationary.

1. Press and hold the *[SELECT]* button.
 - ***MACHINE MENU*** is displayed after 8 seconds.
2. Press the *[SELECT]* button repeatedly until '-> Settings' is displayed.
3. Press the *[START/PULSE]* button.
 - 'End beep = on' or 'End beep = off' is displayed.
4. Press the *[SELECT]* button repeatedly until 'Lid AutoOpen=off' or 'Lid AutoOpen=on' is displayed.
5. Use *[Adjustment buttons]* to set 'off' or 'on'.
 - off: Lid does not unlock automatically.
 - on: Lid unlocks automatically.
6. Press the *[START/PULSE]* button.
 - The setting is stored.
 - 'Store setting...' is displayed briefly.
 - '-> Settings' is then displayed.

7.  Press the *[STOP/OPEN]* button once to exit the '-> Settings' menu
or
Press the *[STOP/OPEN]* button twice to exit the '*MACHINE MENU*'.


7.4.8 Centrifugation data displayed after switching on

The centrifugation data of program 1 or the last program used is displayed after switching on.

1.  Press and hold the *[SELECT]* button.
 - '*MACHINE MENU*' is displayed after 8 seconds.
2.  Press the *[SELECT]* button repeatedly until '-> Settings' is displayed.
3.  Press the *[START/PULSE]* button.
 - 'End beep = on' or 'End beep = off' is displayed.
4.  Press the *[SELECT]* button repeatedly until 'Start Pr = First' or 'Start Pr = Last' is displayed.
5.  Use *[Adjustment buttons]* to set 'Last' or 'First'.
Last = last program used
First = program 1
6.  Press the *[START/PULSE]* button.
 - The setting is stored.
'Store setting...' is displayed briefly.
-> Settings' is then displayed.
7.  Press the *[STOP/OPEN]* button once to exit the '-> Settings' menu
or
Press the *[STOP/OPEN]* button twice to exit the '*MACHINE MENU*'.

7.4.9 Indicator backlight

The indicator backlight can be switched off after 2 minutes to save energy. The rotor is stationary.

1.  Press and hold the *[SELECT]* button.
 - '*MACHINE MENU*' is displayed after 8 seconds.
2.  Press the *[SELECT]* button repeatedly until '-> Settings' is displayed.
3.  Press the *[START/PULSE]* button.
 - 'End beep = on' or 'End beep = off' is displayed.
4.  Press the *[SELECT]* button repeatedly until 'Power save=off' or 'Power save=on' is displayed.
5.  Use *[Adjustment buttons]* to set 'off' or 'on'.
off: Backlight is switched off.
on: Backlight is switched on.
6.  Press the *[START/PULSE]* button.
 - The setting is stored.
'Store setting...' is displayed briefly.
-> Settings' is then displayed.

7. Press the [STOP/OPEN] button once to exit the '-> Settings' menu or
 Press the [STOP/OPEN] button twice to exit the "MACHINE MENU".

8 Cleaning and care

8.1 Overview table

Chap.	Task to execute	if required	daily	weekly	Annually	Page
8	Cleaning and care					35
8.3	Cleaning					36
8.3	Cleaning the device		X			36
8.3	Cleaning accessories			X		36
8.4	Disinfection					36
8.4	Disinfecting the device	X				37
8.4	Disinfecting the accessories	X				37
8.5	Maintenance					37
8.5	Greasing the rubber seal of the centrifuging chamber			X		37
8.5	Trunnion greasing			X		37
8.5	Checking the accessories			X		37
8.5	Centrifuging chamber damage inspection				X	38
8.5	Greasing the motor shaft				X	38
8.5	Accessories with a limited service life	X				38
8.5	Replacing centrifuge tubes	X				38

8.2 Instructions for cleaning and disinfection

DANGER

Risk of contamination for the user due to inadequate cleaning or failure to observe the cleaning instructions.

- Observe cleaning instructions.
- Wear personal protective equipment when cleaning the device.
- Observe laboratory regulations (e.g. TRBAs, the German Protection against Infection Act, hygiene plan) for handling biological agents.

- The device and its accessories must not be cleaned in dishwashers.
- Only perform hand cleaning and liquid disinfection.
- The water temperature must not exceed 25 °C.
- To prevent any corrosion due to use of detergents or disinfectants, it is essential to follow the special application instructions provided by the manufacturers of the detergent or disinfectant.

Disinfectant:

- Surface disinfectant (not disinfectant for hands or instruments)
- Ethanol as the sole active substance.
Do not use an ethanol-propanol mixture to disinfect the viewing window in the lid of the device.
- Concentration not less than 30%
- pH: 6 – 8
- Non-corrosive

8.3 Cleaning

Cleaning the device

1. → Open the lid.
2. → Switch off the device and disconnect it from the power supply.
3. → Remove accessories.
4. → Clean the centrifuge housing and the centrifuging chamber with soap or a mild detergent and a damp cloth.
5. → Remove any detergent residues with a damp cloth after using detergents.
6. → The surfaces must be dried immediately after cleaning.
7. → Dry the centrifuging chamber with an absorbent cloth if condensation forms.

Cleaning accessories

1. → Clean the accessories using the detergent and a damp cloth.
2. → Remove any detergent residues with a damp cloth after using detergents.
3. → Dry the accessories immediately after cleaning using a lint-free cloth and oil-free compressed air. Dry all cavities completely using oil-free compressed air.

8.4 Disinfection

Disinfection must always be preceded by cleaning the components concerned.

See → Chapter 8.3 'Cleaning' on page 36

Disinfectant concentration and application time according to the manufacturer's instructions.

Disinfecting the device

CAUTION

Risk of injury due to ingress of water or other liquids.

- Protect the device against external liquids.
- Do not disinfect the device using spray.

1. Open the lid.
2. Switch off the device and disconnect it from the power supply.
3. Remove accessories.
4. Clean the housing and centrifuging chamber using disinfectant.
5. Remove any disinfectant residues with a damp cloth after using disinfectants.
6. The surfaces must be dried immediately after cleaning.

Disinfecting the accessories

1. Disinfect the accessories using the disinfectant.
2. Wet all cavities with bubble-free disinfectant.
3. Remove the disinfectant residues or leave them to dry after using disinfectants.

Autoclaving

The following accessories may be autoclaved at 121 °C / 250 °F (20 min):

- Swing-out rotors
- Aluminium angle rotors
- Metal buckets
- Lid with bioseal
- Inserting

No statement can be made about the resulting degree of sterility.

The lids of the rotors and bucket must be removed before autoclaving.

Autoclaving accelerates the ageing of materials. It may cause changes in colour. After autoclaving, the rotors and accessories are to be visually inspected for damage and any damaged parts are to be replaced immediately.

The sealing ring in question is to be replaced if there are signs of cracking, embrittlement or wear. For lids with non-replaceable sealing rings, the whole lid must be replaced.

8.5 Maintenance

Greasing the rubber seal of the centrifuging chamber

1. Rub the sealing ring lightly with a rubber care product.

Trunnion greasing

1. Remove accessories.
2. Clean the trunnions.
3. Remove any detergent residues with a damp cloth after using detergents.
4. Grease the trunnions and suspension with Hettich Tubenfett 4051.
5. Excess grease in the centrifuging chamber must be removed.

Checking the accessories

1. The accessories shall be checked for wear and corrosion damage.
2. Check that the rotor is firmly seated.

Centrifuging chamber damage inspection

→ Check the centrifuging chamber for damage.

Greasing the motor shaft

1. → Remove accessories.
2. → Clean the motor shaft.
3. → Remove any detergent residues with a damp cloth after using detergents.
4. → Grease the motor shaft with Hettich Tubenfett 4051.
5. → Excess grease in the centrifuging chamber must be removed.

Accessories with a limited service life

The use of certain accessories is time-limited. For safety reasons, the accessories must no longer be used when either the maximum number of permissible run cycles marked on them or the expiry date marked on them has been reached.

- The maximum permissible number of run cycles or the expiry date can be seen on the accessories.
- The centrifuge is equipped with a cycle counter.

Replacing centrifuge tubes

CAUTION

Risk of injury from broken glass.

Broken glass may cause glass splinters and contaminated liquids to be found inside the centrifuge.

- Wear cut-resistant gloves.
- Wear protective goggles and a face mask.

Broken parts of the tube, glass splinters and spilled centrifuge material must be removed completely in the event of leakage or if a centrifuge tube breaks. Glass splinters that are not removed will cause further glass breakage.

The rubber inserts and the plastic sleeves of the rotors must be replaced after a glass breakage.

Disinfection must be carried out if the material is infectious.

9 Troubleshooting

9.1 Fault description

Customer service must be notified if the fault cannot be rectified based on the fault table. State the centrifuge type and serial number. Both numbers can be seen on the type plate of the centrifuge.

* Error number does not appear on the display.

Fault description	Cause	Remedy
no display	No power. Mains input fuses defective.	<ul style="list-style-type: none"> ■ Check the supply voltage. ■ Check the mains input fuse. ■ The mains switch is in switch position <i>///</i>
IMBALANCE	The rotor is unevenly loaded.	<ul style="list-style-type: none"> ■ Open the lid. ■ Check the loading of the rotor. ■ Repeat the centrifugation run.

Fault description	Cause	Remedy
RPM > ROTOR MAX	Speed in the selected program greater than the maximum rotor speed.	<ul style="list-style-type: none"> Check and correct the speed.
MAINS INTERRUPT	Loss of mains power during the centrifugation run. The centrifugation run was not completed.	<ul style="list-style-type: none"> Open the lid. Press the <i>[START/PULSE]</i> button. If required: Repeat the centrifugation run.
R... WRONG ROTOR (R = rotor code)	The rotor being used is not approved for the device.	<ul style="list-style-type: none"> Install a rotor that is approved for the device.
KEYBOARD-ERROR	Error/defect in electronics.	<ul style="list-style-type: none"> Perform a MAINS RESET.
TACHO - ERROR 1, 2, 96.1	Speed pulse failure. No rotor inserted. Electronics defective.	<ul style="list-style-type: none"> Perform a MAINS RESET after the displayed waiting time (150 seconds) has elapsed. Check whether the rotor is inserted.
LID ERROR 4.1-4.127	Lid lock error.	<ul style="list-style-type: none"> Perform a MAINS RESET.
OVER SPEED 5	Overspeed.	<ul style="list-style-type: none"> Perform a MAINS RESET.
ROTOR ERROR 10.1-10.6	Rotor coding error.	<ul style="list-style-type: none"> Perform a MAINS RESET.
VERSION ERROR 12	Wrong centrifuge model detected. Error / defect in the electronics.	<ul style="list-style-type: none"> Perform a MAINS RESET.
UNDER SPEED 13	Underspeed.	<ul style="list-style-type: none"> Perform a MAINS RESET.
CTRL ERROR 22-25.4	Error/defect in electronics.	<ul style="list-style-type: none"> Perform a MAINS RESET.
CRC ERROR 27, 27.1	Error/defect in electronics.	<ul style="list-style-type: none"> Perform a MAINS RESET.
COM ERROR 31-36	Error/defect in electronics.	<ul style="list-style-type: none"> Perform a MAINS RESET.
FC ERROR 60, 61.1-61.21, 61.64-61.142	Error/defect in electronics.	<ul style="list-style-type: none"> Perform a MAINS RESET.
TACHO ERR 61.22	Speed measurement error.	<ul style="list-style-type: none"> The device must not be switched off while the message "Wait" is displayed. Perform a MAINS RESET when the message "wait ...!" is no longer displayed.
FC ERROR 61.23	Speed measurement error.	<ul style="list-style-type: none"> The device must not be switched off while the message "Wait" is displayed. Perform a MAINS RESET when the message "wait ...!" is no longer displayed.
FC ERROR 61.153	Error/defect in electronics.	<ul style="list-style-type: none"> Perform a MAINS RESET. Check the loading of the rotor. Repeat the centrifugation run.

Fault description	Cause	Remedy
VERS. ERR 61.154	Invalid machine version.	■ Perform a MAINS RESET.

9.2 Performing a MAINS RESET

1. ➔ Set the mains switch to [0].
2. ➔ Wait 10 seconds.
3. ➔ Set the mains switch to [I].

9.3 Emergency release

The lid cannot be unlocked by the motor in the event of a power failure. Emergency unlocking by hand must be performed.

! WARNING

Risk of electric shock due to maintenance and servicing work on live device.

- Disconnect the device from the mains before carrying out repairs and maintenance.

WARNING

Danger of cutting and crushing due to moving rotor.

- Do not open the lid until the rotor has stopped.

Personnel:

- Trained user

1. ➔ Look through the window in the lid to ensure that the rotor is stationary.
2. ➔ Insert the hex key horizontally into the hole (1) and turn anticlockwise until the lid opens.
3. ➔ Remove the hex key from the hole (1).


Fig. 20: Emergency release

1 Hole

9.4 Replacing the mains input fuse

! WARNING

Risk of electric shock due to maintenance and servicing work on live device.

- Disconnect the device from the mains before carrying out repairs and maintenance.

Personnel:

- Trained user


Fig. 21: Mains input fuse

- 1 Fuse holder
- 2 Snap lock

The mains fuses are located next to the mains switch.

The mains switch is in switch position [O]

1. ➤ Disconnect the mains cable from the device plug.
2. ➤ Press the snap locks (2) against the fuse holder (1) and pull them out.
3. ➤ Replace the defective mains input fuses.
Only use fuses with the nominal value specified for the type: see the table below.
4. ➤ Push in the fuse holder (1) until the snap lock engages.
5. ➤ Reconnect the device to the mains.

Model	Type	Fuse	Order no.
EBA 280	1101	T 3.15 AH/250 V	E997
EBA 280	1101-01	T 6.3 AH/250 V	2266
EBA 280 S	1102	T 3.15 AH/250 V	E997
EBA 280 S	1102-01	T 6.3 AH/250 V	2266

10 Disposal

10.1 General instructions

The device can be disposed of via the manufacturer.

A Return Material Authorisation (RMA) form must always be requested for a return.

If necessary, contact the Technical Service Department of the manufacturer:

- **Andreas Hettich GmbH & Co. KG**
- Föhrenstraße 12
- 78532 Tuttlingen, Germany
- Phone: +49 7461 705 1400
- E-Mail: service@hettichlab.com

! WARNING

Risk of pollution and contamination for people and the environment.

When disposing of the centrifuge, people and the environment may be polluted or contaminated by incorrect or improper disposal.

- Removal and disposal may be carried out only by a trained and authorized service personnel.

The device is intended for the commercial sector ("Business to Business" - B2B).

According to Directive 2012/19/EU, the devices may no longer be disposed of with household waste.

The appliances are assigned to the following groups according to the Stiftung Elektro-Altgeräte Register (EAR (German foundation under civil law)):

Fig. 22: Household waste ban

- Group 5 (small appliances)

The crossed-out wheeled bin symbol indicates that the device must not be disposed of with household waste. Regulations governing disposal of such devices may differ in individual countries. If necessary, contact the supplier.

11 Index

A

Accessories.	14
biosafety systems.	36
disinfection.	37
the rotor.	37
with limited service life.	38
Adjustment during centrifugation.	24
Autoclaving.	37

B

Brake level.	25
Buffer, automatic.	28

C

Care Intervals.	35
Centrifugation in continuous operation.	23
with higher substance density.	27
with time preselection.	23
Centrifugation data after switching on.	34
Centrifuging chamber inspection.	38
Centrifuging radius RAD.	26
Cleaning.	36
Cleaning and disinfection Instructions.	35
Connecting the centrifuge.	18
Continuous operation.	23
Cycle counter.	30
Entering a maximum value.	30

D

Device biosafety systems.	36
disinfection.	37
Disinfection.	36
Disposal.	41

E

Enabling/disabling an audible signal.	32
dual time.	31
Error messages.	38

F

Fastening removal.	17
the transport lock.	16
Filling.	21
Foreseeable misuse.	6

G

General safety instructions.	7
Greasing the rubber seal.	37, 38

I

Intended use.	5
-----------------------	---

L

Labels on the device.	12
on the packaging.	11

Lid

the lid.	19
------------------	----

Loading.	21
------------------	----

M

MAINS RESET.	40
Maintenance.	37
Intervals.	35

N

Non-intended use.	6
---------------------------	---

O

Opening programs.	27, 28
Operator's responsibility.	7
Original spare parts.	14

P

Personal protective equipment.	6
Personnel qualifications.	6
Personnel training.	7
Protective equipment.	6

Q

Querying centrifugation runs.	31
operating hours.	31
system information.	29

R

Ramp-up level.	24
Rating plate.	10
Relative centrifugal force RCF.	26
Removing the rotor.	20, 21, 22
Replacing centrifuge tubes.	38
Reset counter Cycle counter.	30
Return.	15
Rotor detection.	29
Runtime programs.	25
Start of timing.	25

S

Safety instructions.	7
Saving programs.	28
Scope of delivery.	14
Setting up the centrifuge.	17
Short-term centrifugation.	24

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

- 1) Sample temp. rise during maximum speed and 1 hour running time
- 2) Observe the tube manufacturer's instructions




1142	1127-A						
Ausschwingrotor 12-fach / Swing out rotor 12-times  $\angle 60^\circ$							
	0553						
							
Kapazität / capacity ml	5	1,6 - 5	2,6 - 3,4	2,7 - 3			
Maße / dimensions $\varnothing \times L$ mm	13 x 75	13 x 75	13 x 65	11 x 66			
Anzahl p. Rotor / number p. rotor	12						
Drehzahl / speed RPM	5000						
RZB / RCF ²⁾	2963						
Radius / radius mm	106						
 9 (97%) sec	10						
 9 sec	12						
Probenerwärmung/Sample temp. rise K ¹⁾	10						

1148	1131-A				1132-A			
Ausschwingrotor 8-fach / Swing out rotor 8-times  $\angle 90^\circ$								
	0553				0501		2079	
								
Kapazität / capacity ml	5	1,6 - 5	2,6 - 3,4	2,7 - 3	6	4 - 5,5	10	4 - 7
Maße / dimensions $\varnothing \times L$ mm	13 x 75	13 x 75	13 x 65	11 x 66	12 x 82	15 x 75	17 x 70	16 x 75
Anzahl p. Rotor / number p. rotor	8							
Drehzahl / speed RPM	5000							
RZB / RCF ²⁾	2991							
Radius / radius mm	107							
 9 (97%) sec	8							
 9 sec	10							
Probenerwärmung/Sample temp. rise K ¹⁾	9							

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




- 1) Sample temp. rise during maximum speed and 1 hour running time
- 2) Observe the tube manufacturer's instructions.
- 3) without lid

1133							
Winkelrotor 12-fach / (mit Dekantierhilfe) Angle rotor 12-times (with decanter aid)  ∠ 35°							
		0553	0501	0578			
							
Kapazität / capacity	ml	5	6	7			
Maße / dimensions	∅ x L mm	12 x 75	12 x 82	12 x 100			
Anzahl p. Rotor / number p. rotor		12					
Drehzahl / speed	RPM	5000					
RZB / RCF	²⁾	2879					
Radius / radius	mm	103					
 g (97%)	sec	8					
 g	sec	10					
Probenerwärmung/Sample temp. rise	K ¹⁾	18					

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- 1) Sample temp. rise during maximum speed and 1 hour running time
- 2) Observe the tube manufacturer's instructions.

LABORATORY EQUIPMENT MAINTENANCE, REPAIR, CALIBRATION AND SALES

Established in 1987, Henderson Biomedical is the UK's leading laboratory equipment sales and service provider. Our knowledgeable team can provide you with excellent sales advice on a range of different types of laboratory equipment including centrifuges, refrigerators, freezers and heat sealers.

Henderson Biomedical is also able to provide you with first class after-sales service and calibration of your laboratory equipment. We are an **ISO 17025 (UKAS) accredited calibration test laboratory** and our team of Field Service Engineers cover the whole of the United Kingdom.

Please contact us for more information on the types of equipment we supply and the different after-sales services we can offer.

Henderson Biomedical

Unit 3, Swan Close
Croydon CR0 2DZ
United Kingdom

Tel: 020 8663 4610

For sales enquiries: sales@henderson-biomedical.co.uk
For all other enquiries: info@henderson-biomedical.co.uk
www.henderson-biomedical.co.uk