

memmert

trust the best



Waterbath

WTB

accessories

**Intuitive.
Flexible.
Calibratable.**

www.memmert.com

Waterbath WTB accessories

IDEAL FOR THE TEMPERATURE CONTROL OF OINTMENTS, EMULSIONS, SAMPLES, PLATES, AND NUTRIENT SOLUTIONS IN THE LABORATORY, AS WELL AS WARM STORAGE AND CALIBRATION IN THE INDUSTRIAL SEGMENT

Flexible and modularly expandable

The waterbath WTB adapts perfectly to your needs with its modular concept. When required by an application, a shaking device, Peltier cooling unit, or a circulation pump can be easily retrofitted.

Keeps the load in motion

If the liquid in the reaction vessels must not rest, the shaking device SV ensures constant movement in the horizontal direction. The shaking frequency, which comes with a range between 30 and 150 strokes per minute, can be set on the display of the shaking device, program runtime and alarm protection.



Active cooling with environmentally friendly Peltier technology

Set the desired temperature directly on the touchscreen of the waterbath WTB. With the Peltier cooling device CDP115, the Memmert waterbath reaches a temperature of up to +10 °C in all sizes. The water is permanently circulated and cooled in the process using special water cooling hoses.

Improved thermal homogeneity

The pump device, CPM650, ensures continuous circulation of the temperature control medium. This self-priming pump is suitable for continuous operation and can be flexibly switched on or off at any time. The flow capacity is 1.2 liters per minute.



Waterbath cover, bottom shelf, and racks

To ensure efficient and energy-saving work, we recommend using a gable or flat cover. Special covers, with corresponding cutouts, as well as various inserts and racks for test tubes and Erlenmeyer flasks, ensure that your load fits perfectly. We also recommend using a bottom shelf so that your load is not placed directly on the floor. With the bottom shelf, the temperature control medium are heated evenly from all sides.





LABORATORY EQUIPMENT MAINTENANCE, REPAIR, CALIBRATION AND SALES

Established in 1987, Henderson Biomedical is the UK's leading laboratory equipment sales and service provider. Our knowledgeable team can provide you with excellent sales advice on a range of different types of laboratory equipment including centrifuges, refrigerators, freezers and heat sealers.

Henderson Biomedical is also able to provide you with first class after-sales service and calibration of your laboratory equipment. We are an **ISO 17025 (UKAS) accredited calibration test laboratory** and our team of Field Service Engineers cover the whole of the United Kingdom.

Please contact us for more information on the types of equipment we supply and the different after-sales services we can offer.

Henderson Biomedical

Unit 3, Swan Close
Croydon CR0 2DZ
United Kingdom

Tel: 020 8663 4610

**For sales enquiries: sales@henderson-biomedical.co.uk
For all other enquiries: info@henderson-biomedical.co.uk
www.henderson-biomedical.co.uk**

