

memmert

trust the best

NEW



Peltier-cooled
incubator

**IPP
eco**

**Efficient.
Homogeneous.
Reliable.**

Peltier-cooled incubator IPPeco

IDEAL FOR CULTIVATION ABOVE AND BELOW ROOM TEMPERATURE,
DETERMINING COLONY COUNT OR DURABILITY TESTS



Maximum energy efficiency compared to cooled incubators with a compressor.



Absolutely homogeneous distribution of temperature.



Greater range of application thanks to cooling down to the freezing point.



Always a good climate: refrigerant-free, smooth running, compact, low-maintenance.

Model sizes: 30, 55, 110, 260, 410, 750, 1060, 1400, 2200
Temperature adjustment ranges: 0 to +70 °C (IPP110eco - IPP2200eco)


Advanced Peltier Technology

As a technology leader for Peltier laboratory equipment, we set new standards: perfectly coordinated systems and excellent efficiency make the Peltier-cooled incubator IPPeco the ultimate when it comes to energy efficiency and reliability.

Homogeneous, constant and reproducible

Despite the significantly reduced energy consumption, the IPPeco's performance is impressive. Heating up, cooling down and recovering after the door is opened occur fine-tuned and yet at top speed. The high-precision control guarantees an absolutely homogeneous and stable temperature distribution.

Appliance comparison shows tremendous potential for energy savings





LABORATORY EQUIPMENT MAINTENANCE, REPAIR, CALIBRATION AND SALES

Established in 1987, Henderson Biomedical is the UK's leading laboratory equipment sales and service provider. Our knowledgeable team can provide you with excellent sales advice on a range of different types of laboratory equipment including centrifuges, refrigerators, freezers and heat sealers.

Henderson Biomedical is also able to provide you with first class after-sales service and calibration of your laboratory equipment. We are an **ISO 17025 (UKAS) accredited calibration test laboratory** and our team of Field Service Engineers cover the whole of the United Kingdom.

Please contact us for more information on the types of equipment we supply and the different after-sales services we can offer.

Henderson Biomedical

Unit 3, Swan Close
Croydon CR0 2DZ
United Kingdom

Tel: 020 8663 4610

For sales enquiries: sales@henderson-biomedical.co.uk
For all other enquiries: info@henderson-biomedical.co.uk
www.henderson-biomedical.co.uk

