

Floorstanding centrifuges

ROTO SILENTA 630 RS

Your first choice up to 12 liters

The ROTO SILENTA 630 RS is highly rated in the rapid and efficient centrifugation of blood bags. It is able to spin up to 12 liters of sample or 12 blood bags per run. Four and six-place rotor options are also available to offer the best possible fit. In industrial or R&D settings the ROTO SILENTA accommodates six vessels up to 2,000 ml each.

This unit includes refrigeration with a temperature range from -20°C to +40°C.



Highlights

- RPM: 50 - 6,000 min⁻¹
Adjustable in increments of 10
- RCF: 50 - 6,520
Adjustable in increments of 1
- max. Capacity: 6 x 2,000 ml
- the first choice in blood banking centrifuges
- choice of 3 rotors
- Medical device according to regulation (EU) 2017/745
- easy operation with keypad and control knob
- 89 program memories for more individuality
- 9 acceleration + 19 deceleration stages
- coolable from -20 bis +40 °C with pre-cooling function
- data reporting system (optional)

Features

- metal housing and lid
- large viewing port in the lid
- powered lid lock
- lid dropping protection
- emergency lid lock release
- stainless steel chamber
- automatic rotor recognition
- brushless drive
- error display
- display in °C and °F possible
- imbalance detection and switch-off

Technical data

ROTO SILENTA 630 RS refrigerated

voltage *)	400 V 3 ~ + N
frequency	50 – 60 Hz
consumption	9,700 VA
emission, immunity	EN/IEC 61326-1, class B
max. capacity	6 x 2,000 ml
max. RPM	6,000 min ⁻¹
max. RCF	6,520
running time	1 – 999 min: 59 s, ∞ continuous run
dimensions (WxDxH)	813 x 1,015 x 973 mm
weight	approx. 355 kg
temperature control, infinitely variable	from -20 to +40 °C
Cat. No.	5005

208 – 220 V +6 / -10 % 3- (+N) + PE,
with internal transformer

consumption

weight

5005-08

9,000 VA

approx. 401 kg

*) Other voltages on request.

GMP versions

		Power supply	Frequency	Cat. No.
ROTO SILENTA 630 RS	GMP, internal cooling unit with water-cooled condenser	400 V 3 ~ +N	50–60 Hz	5005-80
ROTO SILENTA 630 RS	GMP, external cooling unit	400 V 3 ~ +N	50–60 Hz	5005-90

Available rotors

SWING-OUT ROTORS

	angle	max. RPM	max. capacity	Cat. No.	page
 swing-out rotor, 6-place	90°	4,500 min ⁻¹	6 x 2,000 ml	4176	3
 swing-out rotor, 4-place	90°	4,500 min ⁻¹	4 x 1,600 ml	4174	3

ANGLE ROTORS

 angle rotor, 6-place	25°	6,000 min ⁻¹	6 x 250 ml	4570	9
--	-----	-------------------------	------------	-------------	---

— Swing-out rotor, 6- / 4-place | 4176 / 4174

					
Rotor			+	Bucket	
max. RPM		4,500 min ⁻¹		Cat. No.	4524-A
max. RCF	Rotor 4176 / 4174	6,498 / 5,683			
max. capacity		12 / 8 blood bag			
run-up run-down, braked in sec		125 197			
angle temperature in °C ¹⁾		90° +16 / +10			
Cat. No.		4176 / 4174			



			
Blood bag			
capacity in ml	500	500	750
blood bag per system	4-place	4-place	1-place
max. RCF ²⁶⁾	Rotor 4176 / 4174	6,498 / 5,683	6,498 / 5,683
radius in mm	287 / 251	287 / 251	287 / 251
			
Insert			
boring Ø x L in mm	-	-	-
blood bag systems per rotor	12 / 8	12 / 8	12 / 8
Cat. No.	4529-AO,-AM,-AU	4592-B	4592-B

+ Additional blood bank accessories can be found on [page 10](#)

— Swing-out rotor, 6- / 4-place | 4176 / 4174

					
Rotor			+	Bucket	
max. RPM		4,500 min ⁻¹		Cat. No.	4546-A
max. RCF	Rotor 4176 / 4174	6,271 / 5,479			
max. capacity		12 / 8 blood bag			
run-up run-down, braked in sec		125 197			
angle temperature in °C ¹⁾		90° +10 / 0			
Cat. No.		4176 / 4174			



	
Blood bag	
capacity in ml	450
blood bag per system	4-place
max. RCF	Rotor 4176 / 4174
radius in mm	277 / 242


	
Insert	
boring Ø x L in mm	-
blood bag systems per rotor	12 / 8
Cat. No.	4559-A

+ Additional blood bank accessories can be found on [page 10](#)

- 1) For cooled versions: Lowest temperature achievable with precooling and max. speed.
- 26) Includes blood bag hanging to prevent erythrocyte accumulation at low speeds. Different hanging heights allow customer-specific adjustments. (4529-AO top, 4529-AM center, 4529-AU bottom). Hanging blood bags may be centrifuged with a maximum RCF of 1,000.

Swing-out rotor, 6- / 4-place | 4176 / 4174

			
Rotor		Bucket	
max. RPM	4,500 min ⁻¹	Cat. No.	4591-A
max. RCF	Rotor 4176 / 4174 6,498 / 5,705		
max. capacity	12 / 8 blood bag		
run-up run-down, braked in sec	125 197		
angle temperature in °C ¹⁾	90° +16 / +10		
Cat. No.	4176 / 4174		



			
Blood bags			
capacity in ml	450	500	750
blood bag per system	3-place	4-place	1-place
max. RCF	Rotor 4176 / 4174 6,271 / 5,705	6,498 / 5,705	6,498 / 5,705
radius in mm	287 / 252	287 / 251	287 / 251
			
Insert			
boring Ø x L in mm	-	-	-
blood bag systems per rotor	12 / 8	12 / 8	12 / 8
Cat. No.	4598-A	4592-B	4592-B

+ Additional blood bank accessories can be found on [page 10](#)

Swing-out rotor, 6- / 4-place | 4176 / 4174

			
Rotor		Bucket	
max. RPM	3,500 min ⁻¹	Cat. No.	4595-C ³¹⁾
max. RCF	Rotor 4176 / 4174 3,848 / 3,328		
max. capacity	6 x 2,000 ml		
angle temperature in °C ¹⁾	90° +16 / -12		
Cat. No.	4176 / 4174		



	
Vessels	
capacity in ml	max. 1,600 ³⁵⁾ 2,000
W x D x H in mm	- 150 x 100 x 180
max. RCF	Rotor 4176 / 4174 3,821 / 3,328 3,848 / -
radius in mm	279 / 243 281 / -
run-up, in sec	95 / 131 95 / -
run-down, braked in sec	131 / 131 41 / -
	
Insert	
W x D x H in mm	145 x 91 x 146 -
vessels per rotor	6 / 4 6 / -
Cat. No.	4596-A -

- 1) For cooled versions: Lowest temperature achievable with precooling and max. speed.
- 24) At temperatures of over +40 °C and / or when not filled to capacity, bottles may warp during centrifugation.
- 31) Adapter for accommodating sample tubes and blood collection tubes in carrier 4595-C on request.
- 35) Depending on the sample.

— Swing-out rotor, 6- / 4-place | 4176 / 4174


			
Rotor		Bucket	
max. RPM	4,500 min ⁻¹	Cat. No.	4523-A
max. RCF Rotor 4176 / 4174	6,520 / 5,705		
max. capacity	6 / 4 blood bag		
run-up run-down, braked in sec	125 197		
angle temperature in °C ¹⁾	90° +9 / +3		
Cat. No.	4176 / 4174		




	-	-	-
			
Blood bags			
capacity in ml	500	750	1,000
blood bag per system	4-place	1-place	1-place
max. RCF Rotor 4176 / 4174	6,475 / 5,660	6,475 / 5,660	6,520 / 5,705
radius in mm	285 / 250	285 / 250	288 / 252
			
Insert			
boring Ø x L in mm	-	-	-
blood bag systems per rotor	6 / 4	6 / 4	6 / 4
Cat. No.	4516-A	4516-A	-

+ Additional blood bank accessories can be found on [page 10](#)

— Swing-out rotor, 6- / 4-place | 4176 / 4174

			
Rotor		Bucket	
max. RPM	4,500 min ⁻¹	lid	5621
max. RCF Rotor 4176 / 4174	6,316 / 5,524	Cat. No.	4547-B
max. capacity	6 / 4 blood bag		
run-up run-down, braked in sec	125 197		
angle temperature in °C ¹⁾	90° +9 / -9		
Cat. No.	4176 / 4174		



	-	-
		
Blood bags		
capacity in ml	500	450
blood bag per system	4-place	3-place
max. RCF Rotor 4176 / 4174	6,316 / 5,524	6,316 / 5,524
radius in mm	279 / 244	279 / 244
		
Insert		
boring Ø x L in mm	-	-
blood bag systems per rotor	6 / 4	6 / 4
Cat. No.	4548	4548

+ Additional blood bank accessories can be found on [page 10](#)

1) For cooled versions: Lowest temperature achievable with precooling and max. speed.

Swing-out rotor, 6- / 4-place | 4176 / 4174

Rotor

max. RPM	4,500 min ⁻¹
max. RCF	Rotor 4176 / 4174 6,294 / 5,501
max. capacity	6/4 x 1,000 ml
run-up run-down, braked in sec	125 197
angle temperature in °C ¹⁾	90° +4 / -11
Cat. No.	4176 / 4174

Bucket

adapter including lid	4255
Cat. No.	4579-A



Vessels

		5	7	9	15	25	50	100	2.6–3.4	2.7–3	4–5.5	4.5–5	4.9	9–10	10	1.6–5
capacity in ml		12x75	12x100	14x100	17x100	24x100	34x100	44x100	13 x 65	11x66	15x75	11x92	13x90	16x92	15x102	13x75
Ø x L in mm		12x75	12x100	14x100	17x100	24x100	34x100	44x100	13 x 65	11x66	15x75	11x92	13x90	16x92	15x102	13x75
max. RCF ²⁾	Rotor 4176 / 4174	5 750 / 4 935	5 750 / 4 935	5 886 / 5 094	5 886 / 5 094	5 615 / 4 845	5 705 / 4 890	5 683 / 4 867	5 750 / 4 935	5 750 / 4 935	5 886 / 5 094	5 750 / 4 935	5 750 / 4 935	5 886 / 5 094	5 886 / 5 094	5 750 / 4 935
radius in mm		254/218	254/218	260/225	260/225	248/214	252/216	251/215	254/218	254/218	260/225	254/218	254/218	260/225	260/225	254/218



Adapter

boring Ø x L in mm	13 x 58	13 x 58	17,5 x 60	17,5 x 60	26 x 72	36 x 79	45 x 78	13,5 x 58	13 x 58	17,5 x 60	13 x 58	13,5 x 58	17,5 x 60	17,5 x 60	13,5 x 58
vessels per rotor	180 / 120	180 / 120	114 / 76	114 / 76	42 / 28	24 / 16	12 / 8	126 / 84	180 / 120	114 / 76	180 / 120	126 / 84	114 / 76	114 / 76	126 / 84
Cat. No.	4433	4433	4434	4434	4438	4439	4442	4435	4433	4434	4433	4435	4434	4434	4435

Vessels

	blood collection / urine tubes			tubes with screw cap						5127 ²⁴⁾	- ²⁴⁾	0554 ²⁴⁾	0512 ²⁴⁾	4239 ²⁴⁾	4255	Falcon ^{3.15)}
capacity in ml	4-7	4-7	8.5-10	15	50	25	30	10	250	290	650	750	1,000	1,000	1,000	175
Ø x L in mm	13x100	16x75	16x100	17x120	29x115	25x90	25 x 110	16 x 80	61 x 122	62 x 137	97 x 139	97 x 152	96 x 176	98 x 138	98 x 138	61 x 118
max. RCF ²⁾	5,750 / 4,935	5,886 / 5,094	5,886 / 5,094	6,022 / 5,207	5,999 / 5,207	5,818 / 5,026	5,818 / 5,026	5,886 / 5,094	6,113 / 5,320	6,113 / 5,320	5,999 / 5,184	6,294 / 5,501	6,294 / 5,501	6,294 / 5,501	6,294 / 5,501	6,294 / 5,501
radius in mm	254 / 218	260 / 225	260 / 225	266 / 230	265 / 230	257 / 222	257 / 222	260 / 225	270 / 235	270 / 235	265 / 229	278 / 243	278 / 243	278 / 243	278 / 243	278 / 243



Adapter

boring Ø x L in mm	13,5 x 58	17,5 x 60	17,5 x 60	17 x 86	30 x 87	26 x 72	26 x 72	17,5 x 60	62 x 92	62 x 92	98 x 138	98 x 138	98 x 138	-	-	61 x 118
vessels per rotor	126 / 84	114 / 76	114 / 76	72 / 48	30 / 20	42 / 28	42 / 28	114 / 76	6 / 4	6 / 4	6 / 4	6 / 4	6 / 4	6 / 4	6 / 4	6 / 4
Cat. No.	4435	4434	4434	4437	4441	4438	4438	4434	4443	4443	4258	-	-	-	-	4440

Vessels

	Nalgene ^{3.15)}	Nuno ^{3.15)}	Falcon ^{3.15)}	Corning ^{3.15)}	Corning ^{3.15)}
capacity in ml	175	200	225	250	500
Ø x L in mm	61.5x144	60 x 130	61 x 137	60 x 162	96 x 147
max. RCF ²⁾	6,294 / 5,501	6,294 / 5,501	6,294 / 5,501	6,294 / 5,501	6,294 / 5,501
radius in mm	278 / 243	278 / 243	278 / 243	278 / 243	278 / 243



Adapter

boring Ø x L in mm	-	-	-	-	-
vessels per rotor	6 / 4	6 / 4	6 / 4	6 / 4	6 / 4
Cat. No.	4430	4430	4440	4430	4449

- For cooled versions: Lowest temperature achievable with precooling and max. speed.
- Please note that the RCF values indicated refer only to rotor performance. The max. permissible RCF of tubes used should be verified with the individual manufacturers.
The max. RCF for glass tubes annotated with footnote 2) is 4,000.
- When using these tubes, carrier 4255 cannot be closed with its lid.
- At temperatures of over +40 °C and / or when not filled to capacity, bottles may warp during centrifugation.

Swing-out rotor, 6- / 4-place | 4176 / 4174



Rotor

max. RPM	4,500 min ⁻¹
max. RCF	Rotor 4176 / 4174 5,999 / 5,184
max. capacity	6 x 750 ml
run-up run-down, braked in sec	125 197
angle temperature in °C ¹⁾	90° +14 / -1


Cat. No. 4176 / 4174

Bucket

Cat. No.



4522-A



Vessels

	4	5	6	7	12	15	25	50	100	2.6-3.4	2.7-3	4-5.5	4.5-5	7.5-8.2	9-10
capacity in ml	4	5	6	7	12	15	25	50	100	2.6-3.4	2.7-3	4-5.5	4.5-5	7.5-8.2	9-10
Ø x L in mm	10 x 88	12 x 75	12 x 82	12 x 100	16 x 101	17 x 100	24 x 100	34 x 100	40 x 115	13 x 65	11 x 66	15 x 75	11 x 92	15 x 92	16 x 92
max. RCF	5,818 / 5,003	5,818 / 5,003	5,818 / 5,003	5,818 / 5,003	5,818 / 5,003	5,818 / 5,003	5,818 / 5,003	5,818 / 5,003	5,818 / 5,003	5,818 / 5,003	5,818 / 5,003	5,818 / 5,003	5,818 / 5,003	5,818 / 5,003	5,818 / 5,003
radius in mm	257 / 221	257 / 221	257 / 221	257 / 221	257 / 221	257 / 221	257 / 221	257 / 221	257 / 221	257 / 221	257 / 221	257 / 221	257 / 221	257 / 221	257 / 221



Adapter

	11 x 70	12.5 x 32	12.5 x 70	12.5 x 70	16 x 70	17.5 x 70	26 x 70	35 x 70	41.5 x 70	13.2 x 32	12.5 x 32	17.5 x 32	12.5 x 70	17.5 x 70	17.5 x 70
boring Ø x L in mm	11 x 70	12.5 x 32	12.5 x 70	12.5 x 70	16 x 70	17.5 x 70	26 x 70	35 x 70	41.5 x 70	13.2 x 32	12.5 x 32	17.5 x 32	12.5 x 70	17.5 x 70	17.5 x 70
vessels per rotor	378 / 252	288 / 192	288 / 192	288 / 192	150 / 100	180 / 120	66 / 44	36 / 24	24 / 16	180 / 120	288 / 192	180 / 120	288 / 192	180 / 120	96 / 64
Cat. No.	4224	4213-93	4213	4213	4223	4214	4215	4216	4218	4222-93	4213-93	4214-93	4213	4214	4220



Adapter

	17.5 x 70	13.2 x 32	17.5 x 32	13.2 x 70	16 x 70	17.5 x 70	17 x 70	30 x 70	17.5 x 70	26 x 70	30 x 96	62 x 90	62 x 90	94 x 105	97.5 x 105
boring Ø x L in mm	17.5 x 70	13.2 x 32	17.5 x 32	13.2 x 70	16 x 70	17.5 x 70	17 x 70	30 x 70	17.5 x 70	26 x 70	30 x 96	62 x 90	62 x 90	94 x 105	97.5 x 105
vessels per rotor	180 / 120	180 / 120	180 / 120	180 / 120	150 / 100	180 / 120	138 / 92	48 / 32	96 / 64	66 / 44	36 / 24	6 / 4	6 / 4	6 / 4	6 / 4
Cat. No.	4214	4222-93	4214-93	4222	4223	4214	4232	4245-A	4220	4215	4249	4238	4238	4233	4258



Vessels

	750	750	250	500
capacity in ml	750	750	250	500
Ø x L in mm	97 x 152	96 x 135	60 x 162	96 x 147
max. RCF	5,999 / 5,184	5,999 / 5,184	5,818 / 5,003	5,999 / 5,184
radius in mm	265 / 229	265 / 229	257 / 221	265 / 229





Adapter

	97.5 x 105	97.5 x 105	60 x 125	97.5 x 105
boring Ø x L in mm	97.5 x 105	97.5 x 105	60 x 125	97.5 x 105
vessels per rotor	6 / 4	6 / 4	6 / 4	6 / 4
Cat. No.	4258	4258	6322	4258

- 1) For cooled versions: Lowest temperature achievable with precooling and max. speed.
- 2) Please note that the RCF values indicated refer only to rotor performance. The max. permissible RCF of tubes used should be verified with the individual manufacturers. The max. RCF for glass tubes annotated with footnote 2) is 4,000.
- 5) Tested by the TÜV in conformity with DIN EN 61010, section 2 - 020.
- 24) At temperatures of over +40 °C and / or when not filled to capacity, bottles may warp during centrifugation.

Swing-out rotor, 6- / 4-place | 4176 / 4174



Rotor

max. RPM	4,500 min ⁻¹
max. RCF Rotor 4176 / 4174	6,384 / 5,592
max. capacity	6 x 1,000 ml
run-up run-down, braked in sec	125 197
angle temperature in °C ¹⁾	90° +9 / -9
Cat. No.	4176 / 4174

Bucket

lid	5621
Cat. No.	4547-B

Vessels

		5	7	9	15	25	50	100	2.6–3.4	2.7–3	4–5.5	4.5–5	4.9	7.5–8.2	9–10	10
capacity in ml		12x75	12x100	14x100	17x100	24x100	34x100	44x100	13 x 65	11x66	15x75	11x92	13x90	15x92	16x92	15x102
Ø x L in mm		12x75	12x100	14x100	17x100	24x100	34x100	44x100	13 x 65	11x66	15x75	11x92	13x90	15x92	16x92	15x102
max. RCF ²⁾	Rotor 4176 / 4174	5,841 / 5,026	5,841 / 5,026	5,977 / 5,184	5,977 / 5,184	5,728 / 4,913	5,773 / 4,981	5,750 / 4,958	5,841 / 5,026	5,841 / 5,026	5,977 / 5,184	5,841 / 5,026	5,841 / 5,026	5,977 / 5,184	5,977 / 5,184	5,977 / 5,184
radius in mm		258 / 222	258 / 222	264 / 229	264 / 229	253 / 217	255 / 220	254 / 219	258 / 222	258 / 222	264 / 229	258 / 222	258 / 222	264 / 229	264 / 229	264 / 229



Adapter

boring Ø x L in mm	13 x 58	13 x 58	17,5 x 60	17,5 x 60	26 x 72	36 x 79	45 x 78	13,5 x 58	13 x 58	17,5 x 60	13 x 58	13,5 x 58	17,5 x 60	17,5 x 60	17,5 x 60
vessels per rotor	180 / 120	180 / 120	114 / 76	114 / 76	42 / 28	24 / 16	12 / 8	126 / 84	180 / 120	114 / 76	180 / 120	126 / 84	114 / 76	114 / 76	114 / 76
Cat. No.	4433	4433	4434	4434	4438	4439	4442	4435	4433	4434	4433	4435	4434	4434	4434

Vessels

	blood collection / urine tubes						tubes with screw cap				5127 ²⁴⁾	- ²⁴⁾	0554 ²⁴⁾	0512 ²⁴⁾	4239 ²⁴⁾
capacity in ml	1.6–5	4–7	4–7	8	8.5–10	15	50	25	30	10	250	290	650	750	1,000
Ø x L in mm	13x75	13x100	16x75	16 x 125	16x100	17x120	29x115	25x90	25x110	16 x 80	62 x 122	62 x 137	97 x 139	97 x 152	96 x 176
max. RCF ²⁾ Rotor 4176 / 4174	5,841 / 5,026	5,841 / 5,026	5,977 / 5,184	5,977 / 5,184	5,977 / 5,184	6,090 / 5,298	6,090 / 4,958	5,909 / 5,117	5,909 / 5,117	5,977 / 5,184	6,203 / 5,411	6,203 / 5,411	6,384 / 5,592	6,384 / 5,592	6,384 / 5,592
radius in mm	258 / 222	258 / 222	264 / 229	264 / 229	264 / 229	269 / 234	269 / 219	261 / 226	261 / 226	264 / 229	274 / 239	274 / 239	282 / 247	282 / 247	282 / 247



Adapter

boring Ø x L in mm	13,5 x 58	13,5 x 58	17,5 x 60	17,5 x 60	17,5 x 60	17 x 86	30 x 87	26 x 72	26 x 72	17,5 x 60	62x92	62x92	98 x 141	98 x 141	98 x 141
vessels per rotor	126 / 84	126 / 84	114 / 76	114 / 76	114 / 76	72 / 48	30 / 20	42 / 28	42 / 28	114 / 76	6 / 4	6 / 4	6 / 4	6 / 4	6 / 4
Cat. No.	4435	4435	4434	4434	4434	4437	4441	4438	4438	4434	4443	4443	-	-	-

Vessels

	Falcon®	Nalgene®	Nunc®	Falcon®	Corning®	Corning®
capacity in ml	175	175	200	225	250	500
Ø x L in mm	61 x 118	61,5x144,3	60 x 130	61 x 137	60 x 162	96 x 147
max. RCF ²⁾ Rotor 4176 / 4174	6,384 / 5,592	6,384 / 5,592	6,384 / 5,592	6,384 / 5,592	6,384 / 5,592	6,384 / 5,592
radius in mm	282 / 247	282 / 247	282 / 247	282 / 247	282 / 247	282 / 247



Adapter

boring Ø x L in mm	-	-	-	-	-	-
vessels per rotor	6 / 4	6 / 4	6 / 4	6 / 4	6 / 4	6 / 4
Cat. No.	4440	4430	4430	4440	4430	4449


- For cooled versions: Lowest temperature achievable with precooling and max. speed.
- Please note that the RCF values indicated refer only to rotor performance. The max. permissible RCF of tubes used should be verified with the individual manufacturers. The max. RCF for glass tubes annotated with footnote 2) is 4,000.
- 17) When using these tubes, bucket 4547-B cannot be closed with lid 5621.
- 24) At temperatures of over +40 °C and / or when not filled to capacity, bottles may warp during centrifugation.


Swing-out rotor, 6- / 4-place | 4176 / 4174

	
Rotor	
max. RPM	4,500 min ⁻¹
max. RCF Rotor 4176 / 4174	5,999 / 5,184
max. capacity	6 x 750 ml
run-up run-down, braked in sec	125 197
angle temperature in °C ¹⁾	90° +14 / -1
Cat. No.	4176 / 4174

	
Bucket	
Cat. No.	4572



				
Vessels				
capacity in ml	4 – 7	4 – 7	4 – 7	8.5 – 10
Ø x L in mm	13 x 75	13 x 100	16 x 75	16 x 100
max. RCF Rotor 4176 / 4174	5,818 / 5,003	5,818 / 5,003	5,522 / 4,718	5,522 / 4,718
radius in mm	257 / 221	257 / 221	244 / 208	244 / 208

				
Adapter				
boring Ø x L in mm	13.2 x 59	13.2 x 59	17.5 x 59	17.5 x 59
vessels per rotor	336 / 224	336 / 224	210 / 140	210 / 140
Cat. No.	4493	4493	4494	4494

Angle rotor, 6-place | 4570

	
Rotor	
max. RPM max. RCF	6,000 min ⁻¹ 5,594
max. capacity	6 x 250 ml
run-up run-down, braked in sec	64 69
angle max. noise level	25° 54 dB(A)
temperature in °C ¹⁾	-16
Cat. No.	4570

	
Lid bioseal[®] and phenol-resistant	
Cat. No.	INCLUSIVE



	tubes ²⁾			tubes with screw cap								5127 ²⁴⁾
Vessels												
capacity in ml	15	25	94	15	50	10	30	50	85	94	250	
Ø x L in mm	17x100	24x100	38x102	17x120	29x115	16x80	26x95	29x107	38 x 106	38 x 102	61.5x122	
max. RCF ²⁾	5,315	5,152	4,910	5,152	4,830	5,353	5,152	4,830	4,910	4,910	5,594	
radius in mm	132	128	122	128	120	133	128	120	122	122	139	
Adapter												
boring Ø x L in mm	17.6 x 83	26x80	38.6x88	17x106	30x100	16.6x70	26x80	29x90	38.6 x 88	38.6 x 88	61.5x109	
vessels per rotor	42	18	6	30	6	48	18	6	6	6	6	
Cat. No.	5646	5642	5644	5637	5638	5641	5642	5643	5644	5644	-	

1) For cooled versions: Lowest temperature achievable with precooling and max. speed.

2) Please note that the RCF values indicated refer only to rotor performance. The max. permissible RCF of tubes used should be verified with the individual manufacturers. The max. RCF for glass tubes annotated with footnote 2) is 4,000.

5) Tested by the TÜV in conformity with DIN EN 61010, section 2 - 020.

24) At temperatures of over +40 °C and / or when not filled to capacity, bottles may warp during centrifugation.

Accessories for centrifuging blood bag systems



Hettliner

Blood bags can be secured in adapters using HettLiners. After centrifuging, they can be removed by their loops without disturbing the delicate separation of blood components.



Set of compensatory weights

If blood bag carriers are not filled to the same weight, then compensatory weights can be used to make up the difference. One set includes:
40 balancing weights, 2 g each;
20 balancing weights, 5 g each;
10 balancing weights, 10 g each;
5 balancing weights, 20 g each;
5 balancing weights, 40 g each



Balancing inserts

In case there are not enough blood bag systems to occupy every carrier of the rotor, empty carriers can be filled with balancing inserts. Taring weights supplied with the inserts may be used for fine balancing.



Loading station

Secure your samples and avoid accidental spillage with the convenient Loading Station for centrifuge inserts. The Hettich Loading Station allows users to confidently hold and manage their samples without the risk of the inserts toppling over.

Cat. No. 4564

Cat. No. 4566

4584-A for insert 4559-A

4587-A for insert 4592-B

4589-A for insert 4516-A

Cat. No. 4509

Certifications / Registrations

Product certification:




Product registration:



Company certifications:



Dimensions





LABORATORY EQUIPMENT MAINTENANCE, REPAIR, CALIBRATION AND SALES

Established in 1987, Henderson Biomedical is the UK's leading laboratory equipment sales and service provider. Our knowledgeable team can provide you with excellent sales advice on a range of different types of laboratory equipment including centrifuges, refrigerators, freezers and heat sealers.

Henderson Biomedical is also able to provide you with first class after-sales service and calibration of your laboratory equipment. We are an **ISO 17025 (UKAS) accredited calibration test laboratory** and our team of Field Service Engineers cover the whole of the United Kingdom.

Please contact us for more information on the types of equipment we supply and the different after-sales services we can offer.

Henderson Biomedical

Unit 3, Swan Close
Croydon CR0 2DZ
United Kingdom

Tel: 020 8663 4610

For sales enquiries: sales@henderson-biomedical.co.uk
For all other enquiries: info@henderson-biomedical.co.uk
www.henderson-biomedical.co.uk

