



## RLCF0720 Sparkfree Chest Freezer

Part Number	RLCF0720
Capacity (Litres)	215
Style	Chest Freezer
Exterior Dimensions (HxWxD)	865 x 725 x 655mm (including Handle)
Interior Dimensions (HxWxD)	697 x 606 x 536mm
Weight	46kg
Energy Consumption	0.71 kWh/24hr
Operating Temperature	-18°C to -25°C
Factory Set Point	-20°C
Interior Construction	Aluminium
Exterior Construction	White Enamel
Lid Style	Solid
Door Lock	✓
Baskets	1
Basket Dimension (HxWxD)	200 x 250 x 500mm
Feet	4 (adjustable)
Mains Indicator Light	✓
High Temperature Alarm	Visual
Defrost drain	✓
Manual Defrost	✓
Eco Friendly Refrigerant	✓

The **Labcold** Sparkfree range is engineered to store flammable materials safely with no sources of ignition inside the chambers. Anything that could cause a spark is located outside of the chamber.

**Warning: Do not store open containers of volatile substances in this freezer.**

Did you know that 25% of chemical accidents are caused by improper storage?

Source: US Chemical Safety and Laboratory Investigation Board

Designed specifically for laboratory safety with a sparkfree interior and factory fitted with a door lock for optimum security. This model comes with a storage basket to assist with organisation and storage. Please ask at time of order for additional baskets. This chest freezer has manual defrost for greater operator control.



***LABORATORY EQUIPMENT MAINTENANCE, REPAIR, CALIBRATION AND SALES***

Established in 1987, Henderson Biomedical is the UK's leading laboratory equipment sales and service provider. Our knowledgeable team can provide you with excellent sales advice on a range of different types of laboratory equipment including centrifuges, refrigerators, freezers and heat sealers.

Henderson Biomedical is also able to provide you with first class after-sales service and calibration of your laboratory equipment. We are an **ISO 17025 (UKAS) accredited calibration test laboratory** and our team of Field Service Engineers cover the whole of the United Kingdom.

Please contact us for more information on the types of equipment we supply and the different after-sales services we can offer.

**Henderson Biomedical**

Unit 3, Swan Close  
Croydon CR0 2DZ  
United Kingdom

**Tel: 020 8663 4610**

**For sales enquiries: [sales@henderson-biomedical.co.uk](mailto:sales@henderson-biomedical.co.uk)  
For all other enquiries: [info@henderson-biomedical.co.uk](mailto:info@henderson-biomedical.co.uk)  
[www.henderson-biomedical.co.uk](http://www.henderson-biomedical.co.uk)**

