



RLCF2120 Sparkfree Chest Freezer

Part Number	RLCF2120
Capacity (Litres)	607
Style	Chest Freezer
Exterior Dimensions (HxWxD)	865 x 1705 x 655mm (Including Handle)
Interior Dimensions (HxWxD)	697 x 1586 x 536mm
Weight	77kg
Energy Consumption	1.52kWh/24hr
Operating Temperature	-18°C to -25°C
Factory Set Point	-20°C
Interior Construction	Aluminium
Exterior Construction	White Enamel
Lid Style	Solid
Door Hinging	Back (viewed from front)
Door Lock	✓
Baskets	3
Basket Dimension (HxWxD)	200 x 250 x 500mm
Feet	4 (adjustable)
Mains Indicator Light	✓
Defrost Drain	✓
Manual Defrost	✓
Eco Friendly Refrigerant	✓

The **Labcold** Sparkfree range is engineered to store flammable materials safely with no sources of ignition inside the chambers. Anything that could cause a spark is located outside of the chamber.

Warning: Do not store open containers of volatile substances in this freezer.

Did you know that 25% of chemical accidents are caused by improper storage?

Source: US Chemical Safety and Laboratory Investigation Board

Designed specifically for laboratory safety with a sparkfree interior and factory fitted with a door lock for optimum security.

This model comes with 3 storage baskets to assist with organisation and storage. Please ask at time of order for additional baskets.

This chest freezer has manual defrost for greater operator control and a defrost drain to assist with getting rid of any defrost liquids quickly plus it comes with a 1 year UK parts and labour warranty.



LABORATORY EQUIPMENT MAINTENANCE, REPAIR, CALIBRATION AND SALES

Established in 1987, Henderson Biomedical is the UK's leading laboratory equipment sales and service provider. Our knowledgeable team can provide you with excellent sales advice on a range of different types of laboratory equipment including centrifuges, refrigerators, freezers and heat sealers.

Henderson Biomedical is also able to provide you with first class after-sales service and calibration of your laboratory equipment. We are an **ISO 17025 (UKAS) accredited calibration test laboratory** and our team of Field Service Engineers cover the whole of the United Kingdom.

Please contact us for more information on the types of equipment we supply and the different after-sales services we can offer.

Henderson Biomedical

Unit 3, Swan Close
Croydon CR0 2DZ
United Kingdom

Tel: 020 8663 4610

**For sales enquiries: sales@henderson-biomedical.co.uk
For all other enquiries: info@henderson-biomedical.co.uk
www.henderson-biomedical.co.uk**

