



LABCOLD™

medical & scientific refrigeration

RLVL03204 Basic Laboratory Freezer



Constructed from enamelled steel, this freezer from the **Labcold** Basic Range is designed to provide basic storage at a lower price while still providing the confidence of the **Labcold** brand and support.

The freezer has two clear plastic drawers, one which is the interior width of the freezer and a half drawer at the bottom. There is a storage compartment that runs the whole depth of the freezer with a spring loaded flap to keep contents safe at the top. The adjustable feet ensures the freezer can be stabilised and manual defrost provides you with greater operator control over your temperature sensitive contents.

To put you in control, the temperature can be adjusted with an easy dial located on the right hand side of the interior, however, Labcold will by default set the set point to -20°C.

Please note, this freezer is not sparkfree and can not have a lock or digital temperature display fitted to it.

This Basic Freezer is ideal for economic cold storage of items that are not volatile and where precise storage temperature is not product critical. As you would expect from Labcold, this freezer comes backed with a full one years parts and labour warranty.*



Part Number	RLVL03204
Capacity (Litres)	61
Style	Underbench
Exterior Dimensions (HxWxD)	849 x 475 x 448mm
Energy Consumption	0.54 kWh/24h
Noise level	39dB(A)
Temperature Control	✓
Operating Temperature	-20°C
Temperature Set Point Range	-18°C to -23°C
Sparkfree	✗
Interior Construction	Moulded Plastic
Exterior Construction	Enamelled Steel
Drawers	2 + 1 flap
Door Hinging	Right (can be reversed at time of ordering)
Door Style	Solid
Feet	4 adjustable
Manual Defrost	✓
GWP Friendly Refrigerant	✓
Ammonia Free	✓
1 year UK Parts and Labour Warranty	✓

* Labour UK only





LABORATORY EQUIPMENT MAINTENANCE, REPAIR, CALIBRATION AND SALES

Established in 1987, Henderson Biomedical is the UK's leading laboratory equipment sales and service provider. Our knowledgeable team can provide you with excellent sales advice on a range of different types of laboratory equipment including centrifuges, refrigerators, freezers and heat sealers.

Henderson Biomedical is also able to provide you with first class after-sales service and calibration of your laboratory equipment. We are an **ISO 17025 (UKAS) accredited calibration test laboratory** and our team of Field Service Engineers cover the whole of the United Kingdom.

Please contact us for more information on the types of equipment we supply and the different after-sales services we can offer.

Henderson Biomedical

Unit 3, Swan Close
Croydon CR0 2DZ
United Kingdom

Tel: 020 8663 4610

For sales enquiries: sales@henderson-biomedical.co.uk
For all other enquiries: info@henderson-biomedical.co.uk
www.henderson-biomedical.co.uk

