



RLVF0217 Bench Top Sparkfree Freezer

Part Number	RLVF0217
Capacity (Litres)	66
Style	Bench Top
Exterior Dimensions (HxWxD)	740 x 450 x 510mm
Interior Dimensions (HxWxD)	590 x 333 x 300mm
Weight	34kg
Energy Consumption	0.19kWh/24hr
Operating Temperature	-18°C to -25°C
Factory Set Point	-20°C
Interior Construction	Powder Coated Steel
Exterior Construction	Powder Coated Steel
Door Style	Self Closing
Door Hinging	Right (can be reversed please ask at time of order)
Shelves	3
Shelf Dimensions (w x d)	2 of 309 x 275mm 1 of 309 x 150mm
Door Lock	✓
Door Open Alarm	✓
Feet	4
Access Port	✓
Temperature Display	✓
High Temperature Alarm	✓
Low Temperature Alarm	✓
Manual Defrost	✓
Eco Friendly Refrigerant	✓
2 Year Warranty (UK only)	✓

The Labcold Sparkfree range is engineered to store flammable materials safely with no sources of ignition inside the chambers. Anything that could cause a spark is located outside of the chamber.

Warning: Do not store open containers of volatile substances in this freezer.



CLAIM YOUR **2** YEAR WARRANTY

For extra security all sparkfree models are fitted with a built in lock and digital temperature display. The display features high and low audio and visual temperature alarms as standard in addition to a door open alarm.





LABORATORY EQUIPMENT MAINTENANCE, REPAIR, CALIBRATION AND SALES

Established in 1987, Henderson Biomedical is the UK's leading laboratory equipment sales and service provider. Our knowledgeable team can provide you with excellent sales advice on a range of different types of laboratory equipment including centrifuges, refrigerators, freezers and heat sealers.

Henderson Biomedical is also able to provide you with first class after-sales service and calibration of your laboratory equipment. We are an **ISO 17025 (UKAS) accredited calibration test laboratory** and our team of Field Service Engineers cover the whole of the United Kingdom.

Please contact us for more information on the types of equipment we supply and the different after-sales services we can offer.

Henderson Biomedical

Unit 3, Swan Close
Croydon CR0 2DZ
United Kingdom

Tel: 020 8663 4610

**For sales enquiries: sales@henderson-biomedical.co.uk
For all other enquiries: info@henderson-biomedical.co.uk
www.henderson-biomedical.co.uk**

