



LABCOLD™

medical & scientific refrigeration



Labcold -40 to -86°C Ultra Low Temperature 500 litre Freezer

Part Number	LULT0500U
Capacity (Litres)	500
Style	Freestanding
Exterior Dimensions (HxWxD)	1990 x 840 x 880mm
Interior Dimensions (HxWxD)	1310 x 600 x 630mm
Weight	320kg
Operating Temperature	-50°C to -86°C
Interior Construction	Stainless Steel
Exterior Construction	Epoxy Painted Steel
Door Style	Solid
Door Hinging	Right
Remote Alarm contacts	✓
High Temperature Alarm	✓
Low Temperature Alarm	✓
Mains Failure Alarm	✓
Door Open Alarm	✓
Low Battery Alarm	✓
Condenser High Temperature alarm	✓
Digital Temperature Display	✓
Key Operated Door Lock	✗
Password Protected Temperature Settings	✓
Battery Back up (Display and Alarms)	✓
Shelves	3
Inner Doors	4
Castors	4
Access Port	✓
Eco Friendly Refrigerant	✓
5 year warranty 3 x years parts and labour 2 x years parts only	✓

Labcold Ultra Low Temperature freezers use the tried and trusted cascade refrigeration and have been designed to occupy less than 1m² of floor space.



Complete with an easy to use digital controller conveniently situated on the front of the freezer, the temperature settings can be password protected and it always displays the set point for extra peace of mind. It also features alarm battery

back-up so should there be a mains power failure the freezer will alarm immediately to give you as much time as possible to deal with the emergency.

Backed by a 5 year parts and 3 year labour warranty, this Ultra Low Temperature freezer from Labcold is designed to be a robust and reliable place to store your samples.

We can also supply racking, why not ask for details at time of order?





LABORATORY EQUIPMENT MAINTENANCE, REPAIR, CALIBRATION AND SALES

Established in 1987, Henderson Biomedical is the UK's leading laboratory equipment sales and service provider. Our knowledgeable team can provide you with excellent sales advice on a range of different types of laboratory equipment including centrifuges, refrigerators, freezers and heat sealers.

Henderson Biomedical is also able to provide you with first class after-sales service and calibration of your laboratory equipment. We are an **ISO 17025 (UKAS) accredited calibration test laboratory** and our team of Field Service Engineers cover the whole of the United Kingdom.

Please contact us for more information on the types of equipment we supply and the different after-sales services we can offer.

Henderson Biomedical

Unit 3, Swan Close
Croydon CR0 2DZ
United Kingdom

Tel: 020 8663 4610

For sales enquiries: sales@henderson-biomedical.co.uk
For all other enquiries: info@henderson-biomedical.co.uk
www.henderson-biomedical.co.uk

