



LABCOLD™

medical & scientific refrigeration




RLFR2004 Basic Laboratory Refrigerator

Part Number	RLFR2004
Capacity (Litres)	543
Style	Freestanding
Exterior Dimensions (HxWxD)	1960 x 750 x 755mm
Interior Dimensions (HxWxD)	1603 x 622 x 545mm
Interior Construction	ABS Plastic
Exterior Construction	Powder Coated Steel
Net Weight	100kg
Energy Consumption	0.49kWh/24h
Set Point	+4°C
Temperature Range	+2°C to +10°C
Door Style	Solid
Self Closing Door	✓
Door Hinging	Right (can be reversed at time of manufacture)
Door Handle	Factory Fitted
Sparkfree	✗
Shelves	5 shelves 650mm x 530mm
Wheels	2 swivel and 2 swivel with brake
Door Lock	✓
Auto Defrost	✓
Fan Circulated	✓
Access Port	Optional - purchase at time of manufacture only
Digital Temperature Display	✓
High Temperature Alarm	✓
Low Temperature Alarm	✓
Door Open Alarm	✓
Low GWP Refrigerant	✓

With a robust steel exterior and easy clean plastic interior the large refrigerator from the Labcold Basic range is designed to provide laboratory storage at a lower price while still providing the confidence of the Labcold brand and support.

Lower cost production methods have been used in order to offer the best price possible without significant sacrifices in the important areas of function and reliability. This means your refrigerator still comes with a digital display, temperature alarms and a door lock as standard.



This roomy 543 litre refrigerator has five full width shelves and swivel wheels to allow for easy manoeuvrability and ease of cleaning. A digital temperature display shows you the temperature inside the cabinet, plus it has automatic defrost for extra convenience.

As you would expect from Labcold, it is backed by a comprehensive 1 year parts and labour warranty.*

* labour warranty UK only.

01256 705 570

www.labcold.com





LABORATORY EQUIPMENT MAINTENANCE, REPAIR, CALIBRATION AND SALES

Established in 1987, Henderson Biomedical is the UK's leading laboratory equipment sales and service provider. Our knowledgeable team can provide you with excellent sales advice on a range of different types of laboratory equipment including centrifuges, refrigerators, freezers and heat sealers.

Henderson Biomedical is also able to provide you with first class after-sales service and calibration of your laboratory equipment. We are an **ISO 17025 (UKAS) accredited calibration test laboratory** and our team of Field Service Engineers cover the whole of the United Kingdom.

Please contact us for more information on the types of equipment we supply and the different after-sales services we can offer.

Henderson Biomedical

Unit 3, Swan Close
Croydon CR0 2DZ
United Kingdom

Tel: 020 8663 4610

For sales enquiries: sales@henderson-biomedical.co.uk
For all other enquiries: info@henderson-biomedical.co.uk
www.henderson-biomedical.co.uk

