



## RLPL04043 Basic Laboratory Refrigerator

Part Number	RLPL04043
Capacity (Litres)	104
Style	Underbench
Exterior Dimensions (HxWxD)	841 x 495 x 500mm
Energy Consumption	0.323 kWh/24h
Operating Temperature	+4°C
Temperature Set Point Range	+2°C to +10°C
Interior Construction	Moulded Plastic
Exterior Construction	Enamelled Steel
Shelves	3 + 1 bin
Sparkfree	✗
Door Hinging	Right
Door Style	Solid
Feet	2 adjustable feet at front
Internal Light	✓
Auto Defrost	✓
Eco Friendly Refrigerant	✓
Ammonia Free	✓

Constructed from enamelled steel, this fridge from the **Labcold** Basic Range is designed to provide basic storage at a lower price while still providing the confidence of the **Labcold** name and support.

Lower cost production methods have been used in order to offer the best price possible without significant sacrifices in the important areas of function and reliability.



Two adjustable glass shelves.



One glass bin.

This 104 litre fridge has two adjustable shelves and one bin. It comes with an internal light and thermostat and has autodefrost for your convenience.



***LABORATORY EQUIPMENT MAINTENANCE, REPAIR, CALIBRATION AND SALES***

Established in 1987, Henderson Biomedical is the UK's leading laboratory equipment sales and service provider. Our knowledgeable team can provide you with excellent sales advice on a range of different types of laboratory equipment including centrifuges, refrigerators, freezers and heat sealers.

Henderson Biomedical is also able to provide you with first class after-sales service and calibration of your laboratory equipment. We are an **ISO 17025 (UKAS) accredited calibration test laboratory** and our team of Field Service Engineers cover the whole of the United Kingdom.

Please contact us for more information on the types of equipment we supply and the different after-sales services we can offer.

**Henderson Biomedical**

Unit 3, Swan Close  
Croydon CR0 2DZ  
United Kingdom

**Tel: 020 8663 4610**

**For sales enquiries: [sales@henderson-biomedical.co.uk](mailto:sales@henderson-biomedical.co.uk)  
For all other enquiries: [info@henderson-biomedical.co.uk](mailto:info@henderson-biomedical.co.uk)  
[www.henderson-biomedical.co.uk](http://www.henderson-biomedical.co.uk)**

