



LABCOLD™

medical & scientific refrigeration



RLPL04044 Basic Laboratory Refrigerator

Part Number	RLPL04044
Capacity (Litres)	91
Style	Underbench
Exterior Dimensions (HxWxD)	842 x 475 x 448mm
Weight	22kg
Energy Consumption	0.29 kWh/24h
Noise Level	39 dB(A)
Operating Temperature	+4°C
Temperature Set Point Range	+2°C to +10°C
Interior Construction	Moulded Plastic
Exterior Construction	Enamelled Steel
Shelves	2 + 1 bin
Sparkfree	✗
Door Hinging	Right (can be reversed at time of ordering)
Door Style	Solid
Feet	2 adjustable feet at front
Internal Light	✓
Auto Defrost	✓
GWP Friendly Refrigerant	✓
Ammonia Free	✓
1 year UK Parts and Labour Warranty	✓

Constructed from enamelled steel, this refrigerator from the **Labcold** Basic Range is designed to provide basic storage at a lower price while still providing the confidence of the **Labcold** brand and support.

The fridge has two non adjustable shelves and a polycarbonate drawer at the bottom of the refrigerator. The adjustable feet ensures the fridge can be stabilised and with auto defrost and an internal light, the refrigerator has some useful convenience features.

To put you in control, the temperature can be easily adjusted with a dial located on the right hand side of the interior. Labcold will by default set the set point to +4°C.

Please note, this refrigerator is not sparkfree and can not have a lock or digital temperature display fitted to it.

This Basic Refrigerator is ideal for economic cold storage of items that are not volatile and where precise storage temperature is not product critical. As you would expect from Labcold, this Basic Refrigerator comes complete with a full one years parts and labour warranty.*



* Labour UK only



LABORATORY EQUIPMENT MAINTENANCE, REPAIR, CALIBRATION AND SALES

Established in 1987, Henderson Biomedical is the UK's leading laboratory equipment sales and service provider. Our knowledgeable team can provide you with excellent sales advice on a range of different types of laboratory equipment including centrifuges, refrigerators, freezers and heat sealers.

Henderson Biomedical is also able to provide you with first class after-sales service and calibration of your laboratory equipment. We are an **ISO 17025 (UKAS) accredited calibration test laboratory** and our team of Field Service Engineers cover the whole of the United Kingdom.

Please contact us for more information on the types of equipment we supply and the different after-sales services we can offer.

Henderson Biomedical

Unit 3, Swan Close
Croydon CR0 2DZ
United Kingdom

Tel: 020 8663 4610

For sales enquiries: sales@henderson-biomedical.co.uk
For all other enquiries: info@henderson-biomedical.co.uk
www.henderson-biomedical.co.uk

