




RAFR05204 Advanced Under Counter Freezer

Part Number	RAFR05204
Capacity (Litres)	115
Style	Under Counter
Exterior Dimensions (HxWxD)	835 x 600 x 600mm
Energy Consumption	1.5kWh/24h
Operating Temperature	-20°C
Temperature Set Point Range	-18°C to -22°C
Interior Construction	Stainless Steel
Exterior Construction	Stainless Steel
Shelves	2
Shelf dimensions (W x D)	460 x 4075mm 460 x 240mm
Door Hinging	Right (can be reversed at time of order)
Self Closing Door	✓
Door Lock	✓
Levelling Feet	4
Access Port	✓
Remote Alarm Contacts	✓
Digital Temperature Control	✓
Digital Temperature Display	✓
High Temperature alarm	✓
Low Temperature alarm	✓
Door Open alarm	✓
Manual defrost	✓
GWP Friendly Refrigerant	✓
Ammonia Free	✓
2 Year Warranty*	✓

Robustly constructed with a stainless Steel exterior this freezer from the **Labcold** Advanced range is designed to provide close temperature control for temperature sensitive products .



Insulated with Polyurethane Water Base Insulation the refrigerant is also ammonia, HCFC and CFC free and uses a GWP friendly refrigerant.

Fitted with a sturdy door lock, this Advanced freezer is ideal for demanding environments such as busy laboratories and features a manual defrost putting you in control.

Remote alarm contacts and an access port are also fitted as standard.

*Parts only warranty outside of the UK



LABORATORY EQUIPMENT MAINTENANCE, REPAIR, CALIBRATION AND SALES

Established in 1987, Henderson Biomedical is the UK's leading laboratory equipment sales and service provider. Our knowledgeable team can provide you with excellent sales advice on a range of different types of laboratory equipment including centrifuges, refrigerators, freezers and heat sealers.

Henderson Biomedical is also able to provide you with first class after-sales service and calibration of your laboratory equipment. We are an **ISO 17025 (UKAS) accredited calibration test laboratory** and our team of Field Service Engineers cover the whole of the United Kingdom.

Please contact us for more information on the types of equipment we supply and the different after-sales services we can offer.

Henderson Biomedical

Unit 3, Swan Close
Croydon CR0 2DZ
United Kingdom

Tel: 020 8663 4610

**For sales enquiries: sales@henderson-biomedical.co.uk
For all other enquiries: info@henderson-biomedical.co.uk
www.henderson-biomedical.co.uk**

