

hawo

Seal**LINK** | Smart Networking System

SIMPLE INTELLIGENT PACKAGING



SMART NETWORKING SYSTEM

The **innovative** Seal**LINK** networking system takes the packaging process for medical devices to a new, even more efficient level.

The direct interaction of a hawo sealing machine equipped with SealLINK software and our hm 8000 AS-V pouch production machine allows automatic and customised pouch production with just a single scan. Thanks to the intelligent SealLINK technology, pouches that are too short or too long are a thing of the past. The process is continuously optimised with each new instrument scanned by the data matrix reader.



Seal**LINK**
READY

The following models are compatible with SealLINK:


hm 950 DC-V NanoPak
hm 4000 DC-VI
hm 5000 DC-VI PrintPak




Scan the instrument code with the data matrix reader.



SealLINK automatically determines the appropriate pouch length and stores it.



hm 8000 AS-V produces the right pouch directly with the next scan.



Insert the instrument and seal. Done!



Customised pouches up to max. 2,500 mm length.

Seal**LINK** combines work efficiency, process optimisation and material savings in one single system.

Scan one time, store data one time, get the right pouch size each and every time.



LABORATORY EQUIPMENT MAINTENANCE, REPAIR, CALIBRATION AND SALES

Established in 1987, Henderson Biomedical is the UK's leading laboratory equipment sales and service provider. Our knowledgeable team can provide you with excellent sales advice on a range of different types of laboratory equipment including centrifuges, refrigerators, freezers and heat sealers.

Henderson Biomedical is also able to provide you with first class after-sales service and calibration of your laboratory equipment. We are an **ISO 17025 (UKAS) accredited calibration test laboratory** and our team of Field Service Engineers cover the whole of the United Kingdom.

Please contact us for more information on the types of equipment we supply and the different after-sales services we can offer.

Henderson Biomedical

Unit 3, Swan Close
Croydon CR0 2DZ
United Kingdom

Tel: 020 8663 4610

For sales enquiries: sales@henderson-biomedical.co.uk
For all other enquiries: info@henderson-biomedical.co.uk
www.henderson-biomedical.co.uk

