



LABCOLD™


medical & scientific refrigeration



Labcold -80°C Ultra Low Temperature Freezer LULT3585

Using tried and tested cascade refrigeration this quiet and compact Ultra Low Temperature freezer from Labcold features VIP insulation to improve energy usage and maximise internal space without increasing external size.

The controller displays the temperature and features high and low temperature and a door alarm plus it has battery back up so will function in the event of a mains failure.



Intuitive Control Panel

This freezer also has a door lock plus has an access port and remote alarm contacts and can fit up to 5 stainless steel drawers (optional extra) to ensure the contents can be organised.



For extra peace of mind, this freezer has a FREE 5 year parts warranty and a 3 year labour warranty.

Part Number	LULT3585
Capacity (Litres)	35
Style	Desk Top
Exterior Dimensions (HxWxD)	540 x 580 x 765 mm
Interior Dimensions (HxWxD)	280 x 425 x 300 mm
Weight	81kg
Operating Temperature	-80°C
Interior Construction	Stainless Steel
Exterior Construction	Stainless Steel
VIP insulation	✓
Door Style	Solid
Door Hinging	Right
Open Door alarm	✓
High Temperature Alarm	✓
Mains Failure Alarm	✓
Digital Temperature Display	✓
Access Port	✓
Remote Alarm Contacts	✓
Key Operated Door Lock	✓
Battery Back up (Display & Alarms only)	✓
Drawers	5 (optional extra)
Inner Door	✗
Cascade Refrigeration	✓
FREE 5 year parts and 3 year labour warranty (UK only)	✓
Energy Consumption (kWh/24h)	4
Low GWP Refrigerant	✓

01256 705 570

www.labcold.com





LABORATORY EQUIPMENT MAINTENANCE, REPAIR, CALIBRATION AND SALES

Established in 1987, Henderson Biomedical is the UK's leading laboratory equipment sales and service provider. Our knowledgeable team can provide you with excellent sales advice on a range of different types of laboratory equipment including centrifuges, refrigerators, freezers and heat sealers.

Henderson Biomedical is also able to provide you with first class after-sales service and calibration of your laboratory equipment. We are an **ISO 17025 (UKAS) accredited calibration test laboratory** and our team of Field Service Engineers cover the whole of the United Kingdom.

Please contact us for more information on the types of equipment we supply and the different after-sales services we can offer.

Henderson Biomedical

Unit 3, Swan Close
Croydon CR0 2DZ
United Kingdom

Tel: 020 8663 4610

For sales enquiries: sales@henderson-biomedical.co.uk
For all other enquiries: info@henderson-biomedical.co.uk
www.henderson-biomedical.co.uk

