



LABCOLD™

medical & scientific refrigeration



Labcold -80°C Ultra Low Temperature Freezer LULT80100

Using tried and tested cascade refrigeration this new under bench ULT freezer from Labcold also features VIP insulation to keep the footprint even smaller.



Intuitive Control Panel

For sample security the LULT80100 has a door lock plus features all the relevant alarms. It also comes complete with factory fitted remote alarm contacts and an access port.

For extra peace of mind, this freezer has a FREE 5 year parts warranty and a 3 year labour warranty.



**Need racking?
Why not ask for a
quote?**

Part Number	LULT80100
Capacity (Litres)	100
Style	Under bench
Exterior Dimensions (HxWxD)	800 x 950 x 735 mm
Interior Dimensions (HxWxD)	500 x 450 x 450 mm
Weight	120 kg
Operating Temperature	-80°C
Interior Construction	Stainless Steel
Exterior Construction	Stainless Steel
VIP insulation	✓
Door Style	Solid
Door Hinging	Right
Open Door alarm	✓
High Temperature Alarm	✓
Mains Failure Alarm	✓
Digital Temperature Display	✓
Access Port	✓
Remote Alarm Contacts	✓
Key Operated Door Lock	✓
Battery Back up (Display & Alarms only)	✓
Shelves	1
Inner Door	✗
Rollers	✓
Cascade Refrigeration	✓
FREE 5 year parts and 3 year labour warranty (UK only)	✓
Energy Consumption (kWh/24h)	8.7
Eco Friendly Refrigerant	✓

01256 705 570

www.labcold.com





LABORATORY EQUIPMENT MAINTENANCE, REPAIR, CALIBRATION AND SALES

Established in 1987, Henderson Biomedical is the UK's leading laboratory equipment sales and service provider. Our knowledgeable team can provide you with excellent sales advice on a range of different types of laboratory equipment including centrifuges, refrigerators, freezers and heat sealers.

Henderson Biomedical is also able to provide you with first class after-sales service and calibration of your laboratory equipment. We are an **ISO 17025 (UKAS) accredited calibration test laboratory** and our team of Field Service Engineers cover the whole of the United Kingdom.

Please contact us for more information on the types of equipment we supply and the different after-sales services we can offer.

Henderson Biomedical

Unit 3, Swan Close
Croydon CR0 2DZ
United Kingdom

Tel: 020 8663 4610

**For sales enquiries: sales@henderson-biomedical.co.uk
For all other enquiries: info@henderson-biomedical.co.uk
www.henderson-biomedical.co.uk**

