

memmert

trust the best

Waterbath

WWTB

NEW



**Intuitive.
Flexible.
Calibratable.**

www.memmert.com

Waterbath WTB

IDEAL FOR THE TEMPERATURE CONTROL OF OINTMENTS, EMULSIONS, SAMPLES, PLATES AND NUTRIENT SOLUTIONS IN THE LABORATORY, AS WELL AS WARM STORAGE AND CALIBRATION IN THE INDUSTRIAL SEGMENT



Intuitive to operate,
customisable remote screen
CustomView.



Practical and robust:
quick water drainage,
easy to clean.



Can be flexibly
expanded with customised
accessories.



Can be calibrated for
precise temperature settings.

Model sizes: 6, 11, 15, 24, 35, 50

Working temperature range: +5°C above the ambient temperature up to +100°C

Working temperature range with Peltier cooling unit CDP: +10 °C to 100 °C

Intuitive and user friendly

Memmert is known for excellent operating concepts - the new WTB waterbaths are just more proof of this. The touchscreen impresses with its intuitive and convenient menu navigation. A particular highlight: Thanks to the personalisable remote viewing through the CustomView, you always have an eye on the important information.



Practical and robust

Perfect for everyday use in the laboratory: The high-quality stainless steel container with rounded corners is easy to clean and can be completely drained using the standard drain valve without residue left behind. The compact ergonomic design means the Memmert WTB fits on any laboratory table.



Can be expanded in a flexible and modular way

The modular concept of this modern waterbath guarantees maximum flexibility. If the application calls for it, a circulation pump, cooling unit, shaking device and test tube holder or insertable rack can easily be retrofitted.





LABORATORY EQUIPMENT MAINTENANCE, REPAIR, CALIBRATION AND SALES

Established in 1987, Henderson Biomedical is the UK's leading laboratory equipment sales and service provider. Our knowledgeable team can provide you with excellent sales advice on a range of different types of laboratory equipment including centrifuges, refrigerators, freezers and heat sealers.

Henderson Biomedical is also able to provide you with first class after-sales service and calibration of your laboratory equipment. We are an **ISO 17025 (UKAS) accredited calibration test laboratory** and our team of Field Service Engineers cover the whole of the United Kingdom.

Please contact us for more information on the types of equipment we supply and the different after-sales services we can offer.

Henderson Biomedical

Unit 3, Swan Close
Croydon CR0 2DZ
United Kingdom

Tel: 020 8663 4610

**For sales enquiries: sales@henderson-biomedical.co.uk
For all other enquiries: info@henderson-biomedical.co.uk
www.henderson-biomedical.co.uk**

