



LABCOLD™

medical & scientific refrigeration




RLHE1345 -40°C Sparkfree Chest Freezer

Warning: Do not store open containers of volatile substances in this freezer.

Robust chest freezer featuring a spark free interior for laboratory safety.

Fitted with a lid lock as standard, this freezer also comes complete with a digital temperature display and high and low temperature alarms.



Bespoke racking can be provided if required as can baskets at an extra cost. To increase the life of your freezer and to ensure you get optimum performance, preventative maintenance contracts are also available. Why not ask for details when you order?



Digital Temperature Display

Part Number	RLHE1345
Capacity (Litres)	375
Style	Chest Freezer
Exterior Dimensions (HxWxD)	860 x 1500 x 730mm
Interior Dimensions (HxWxD)	640 x 1295 x 445mm
Weight	82kg
Operating Temperature	-20°C to -45°C
Interior Construction	Aluminium
Exterior Construction	Painted Steel
Defrost	Manual
Sparkfree Interior	✓
Door Open Alarm	X
High Temperature Alarm	✓
Low Temperature Alarm	✓
Digital Temperature Display	✓
Key Operated Lid Lock	✓
Feet	4
Energy Consumption	0.8 kWh/24h
Eco Friendly Refrigerant	✓

01256 705 570

www.labcold.com





LABORATORY EQUIPMENT MAINTENANCE, REPAIR, CALIBRATION AND SALES

Established in 1987, Henderson Biomedical is the UK's leading laboratory equipment sales and service provider. Our knowledgeable team can provide you with excellent sales advice on a range of different types of laboratory equipment including centrifuges, refrigerators, freezers and heat sealers.

Henderson Biomedical is also able to provide you with first class after-sales service and calibration of your laboratory equipment. We are an **ISO 17025 (UKAS) accredited calibration test laboratory** and our team of Field Service Engineers cover the whole of the United Kingdom.

Please contact us for more information on the types of equipment we supply and the different after-sales services we can offer.

Henderson Biomedical

Unit 3, Swan Close
Croydon CR0 2DZ
United Kingdom

Tel: 020 8663 4610

For sales enquiries: sales@henderson-biomedical.co.uk
For all other enquiries: info@henderson-biomedical.co.uk
www.henderson-biomedical.co.uk

